



# Welch Allyn Spot Vital Signs LXi

---



---

## Service Manual

**WelchAllyn®**

Advancing Frontline Care™



# Welch Allyn Spot Vital Signs LXi

---

---

## Service Manual

**WelchAllyn®**

Advancing Frontline Care™

Copyright 2012 Welch Allyn. All rights are reserved. To support the intended use of the product described in this publication, the purchaser of the product is permitted to copy this publication, for internal distribution only, from the media provided by Welch Allyn. No other use, reproduction, or distribution of this publication, or any part of it, is permitted without written permission from Welch Allyn. Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Welch Allyn, Spot Vital Signs, SureBP Technology, and SureTemp are registered trademarks of Welch Allyn.

Braun ThermoScan is a registered trademark of the Braun Company.

LNCS is a trademark of, and SET, LNOP, and Masimo are registered trademarks of, Masimo Corporation. Possession or purchase of a Masimo SpO<sub>2</sub>-equipped device does not convey any express or implied license to use the device with unauthorized sensors or cables which would, alone or in combination with this device, fall within the scope of one or more of the patents relating to this device.

Nellcor and *Oxi-Max* are registered trademarks of Nellcor Puritan Bennett Inc.

Health o meter is a registered trademark of Sunbeam Products, Inc. used under license.

Software in this product is Copyright Welch Allyn or its vendors. All rights are reserved. The software is protected by United States of America copyright laws and international treaty provisions applicable worldwide. Under such laws, the licensee is entitled to use the copy of the software incorporated with this instrument as intended in the operation of the product in which it is embedded. The software may not be copied, decompiled, reverse-engineered, disassembled, or otherwise reduced to human-perceivable form. This is not a sale of the software or any copy of the software; all right, title, and ownership of the software remain with Welch Allyn or its vendors.


For information about any Welch Allyn product, call Welch Allyn Technical Support:

|                      |                  |              |                  |
|----------------------|------------------|--------------|------------------|
| USA                  | +1 800 289 2501  | Australia    | +61 2 9638 3000  |
| Canada               | +1 800 561 8797  | China        | +86 21 6327 9631 |
| European Call Center | +353 46 90 67790 | France       | +33 15 5695 849  |
| Germany              | +49 695 098 5132 | Japan        | +81 3 6383 0852  |
| Latin America        | +1 305 669 9003  | Netherlands  | +31 202 061 360  |
| Singapore            | +65 6419 8100    | South Africa | +27 11 777 7555  |
| United Kingdom       | +44 207 365 6780 | Sweden       | +46 85 853 6551  |

## Reorder #4500-89

Manual Material No. 718440 Ver. A

Welch Allyn  
4341 State Street Road  
Skaneateles Falls, NY  
13153 USA



Regulatory Affairs Representative  
Welch Allyn Limited  
Navan Business Park  
Dublin Road  
Navan, County Meath, Republic of Ireland

Printed in USA

|                                     |           |
|-------------------------------------|-----------|
| <b>1 - Safety Summary</b>           | <b>1</b>  |
| Introduction                        | 1         |
| Warnings and Cautions               | 1         |
| General Warnings                    | 1         |
| Blood Pressure Warnings             | 3         |
| Temperature Warnings                | 4         |
| SpO <sub>2</sub> Warnings           | 5         |
| General Cautions                    | 6         |
| Blood Pressure Cautions             | 7         |
| Temperature Cautions                | 7         |
| SpO <sub>2</sub> Cautions           | 7         |
| Electrostatic Discharge (ESD)       | 8         |
| Symbols                             | 9         |
| Safety Symbols                      | 9         |
| Button Symbols                      | 9         |
| Connection Symbols                  | 9         |
| Agency Symbols                      | 10        |
| <b>2 - Overview</b>                 | <b>11</b> |
| Purpose and Scope                   | 11        |
| Other Applicable Documents          | 11        |
| Miscellaneous Mounting Accessories  | 12        |
| Wall Mount Kit                      | 12        |
| Mobile Stand Kit                    | 14        |
| Radio Mounting Accessory            | 16        |
| Barcode Reader Mounting Accessory   | 17        |
| Warranty                            | 18        |
| Spot Vital Signs LXi                | 18        |
| Accessories                         | 18        |
| Service                             | 19        |
| User Serviceable Parts              | 19        |
| Technical Assistance                | 19        |
| Field Replacement Units             | 19        |
| Service Loaners                     | 19        |
| Service Intervals                   | 19        |
| Spot Vital Signs LXi Configurations | 20        |
| Controls                            | 20        |
| Display Window                      | 22        |
| Connections                         | 24        |
| Braun ThermoScan PRO 4000 Lock      | 25        |
| Blood Pressure Hose and Cuff        | 25        |
| Thermometer                         | 26        |
| SpO <sub>2</sub> Sensor             | 27        |
| Radio                               | 27        |
| Quick Reference Card                | 27        |
| AC Power Transformer                | 28        |
| Printer                             | 28        |
| Battery                             | 29        |

|                                                    |           |
|----------------------------------------------------|-----------|
| <b>3 - Functional Overview</b>                     | <b>31</b> |
| Self Test                                          | 31        |
| Internal Configuration Mode                        | 32        |
| <b>4 - Spot Vital Signs LXi Technical Overview</b> | <b>37</b> |
| System Description                                 | 37        |
| Battery System                                     | 37        |
| Battery Charger                                    | 37        |
| "CPU I/O" Power Supply                             | 37        |
| Standby Mode                                       | 38        |
| CPU Core Power Supply                              | 38        |
| +5V Supply                                         | 38        |
| Clock/Calendar Power                               | 38        |
| Mod F NIBP Power                                   | 38        |
| SpO <sub>2</sub> and Thermometer Power             | 38        |
| LCD Power                                          | 38        |
| External Load (RS-232) Power                       | 39        |
| SureTemp Plus Thermometer Interface                | 39        |
| Braun ThermoScan PRO 4000 Thermometer Interface    | 39        |
| SpO <sub>2</sub> Description                       | 39        |
| Radio Description                                  | 40        |
| Bar Code Scanner Description                       | 40        |
| Interface                                          | 40        |
| Printer Description                                | 41        |
| Interconnect Diagram                               | 42        |
| <b>5 - Maintenance and Service</b>                 | <b>43</b> |
| Cleaning                                           | 43        |
| Spot Vital Signs LXi                               | 43        |
| Blood Pressure Cuff                                | 43        |
| Blood Pressure Hose and Cable                      | 43        |
| SureTemp Plus Thermometer                          | 44        |
| Braun ThermoScan PRO 4000 Thermometer              | 44        |
| SpO <sub>2</sub> Sensors                           | 45        |
| Radio                                              | 45        |
| Printer                                            | 45        |
| Battery Replacement                                | 46        |
| Spot Vital Signs LXi                               | 46        |
| Braun ThermoScan PRO 4000                          | 46        |
| Printer                                            | 47        |
| Spot Vital Signs LXi Firmware Upgrade CD           | 48        |
| Spot Vital Signs LXi Radio Configuration CD        | 48        |
| Product Disposal                                   | 48        |
| <b>6 - Calibration and Verification</b>            | <b>49</b> |
| Terms                                              | 49        |
| Required Test Equipment                            | 50        |
| Pneumatic test assembly setup                      | 50        |
| Power Supply Test Setup                            | 52        |

|                                                                                              |                |
|----------------------------------------------------------------------------------------------|----------------|
| Battery Removal / Power Shutdown . . . . .                                                   | 53             |
| Voltage Calibration . . . . .                                                                | 55             |
| Pressure verification test . . . . .                                                         | 57             |
| Blood pressure calibration . . . . .                                                         | 62             |
| NIBP board Initialization . . . . .                                                          | 62             |
| Functional tests . . . . .                                                                   | 67             |
| Temperature Functional Check . . . . .                                                       | 77             |
| Temperature test setup . . . . .                                                             | 77             |
| Changing the Calibration Set Point Temperature on the 9600 Plus Calibration Tester . . . . . | 78             |
| Welch Allyn SureTemp Plus Thermometers . . . . .                                             | 78             |
| Braun ThermoScan Pro 4000 . . . . .                                                          | 82             |
| Changing the Calibration Set Point Temperature . . . . .                                     | 85             |
| <br><b>7 - Specifications . . . . .</b>                                                      | <br><b>89</b>  |
| Performance . . . . .                                                                        | 89             |
| Blood Pressure Accuracy . . . . .                                                            | 89             |
| Temperature Specifications . . . . .                                                         | 90             |
| SpO <sub>2</sub> Specifications . . . . .                                                    | 91             |
| Masimo Sensor Accuracy Guide . . . . .                                                       | 91             |
| Masimo Patents . . . . .                                                                     | 91             |
| Nellcor Sensor Accuracy Guide . . . . .                                                      | 92             |
| Nellcor Patents . . . . .                                                                    | 92             |
| Mechanical . . . . .                                                                         | 93             |
| Electrical . . . . .                                                                         | 93             |
| Patents . . . . .                                                                            | 93             |
| Batteries . . . . .                                                                          | 93             |
| Environmental . . . . .                                                                      | 93             |
| Radio . . . . .                                                                              | 94             |
| Mechanical . . . . .                                                                         | 94             |
| Configuration Options . . . . .                                                              | 94             |
| Guidance and Manufacturer's Declaration For Spot Vital Signs LXi Radio . . . . .             | 95             |
| <br><b>8 - Troubleshooting . . . . .</b>                                                     | <br><b>99</b>  |
| Error Codes . . . . .                                                                        | 99             |
| Radio Test . . . . .                                                                         | 110            |
| Service Work Checklist . . . . .                                                             | 111            |
| Spot LXi Radio Configuration Tool . . . . .                                                  | 112            |
| <br><b>9 - Supplies and Accessories . . . . .</b>                                            | <br><b>113</b> |
| Blood Pressure . . . . .                                                                     | 113            |
| Temperature . . . . .                                                                        | 114            |
| Pulse Oximetry . . . . .                                                                     | 115            |
| Masimo Accessories . . . . .                                                                 | 115            |
| Nellcor Accessories . . . . .                                                                | 117            |
| Miscellaneous . . . . .                                                                      | 118            |
| Service Contracts . . . . .                                                                  | 119            |

|                                               |            |
|-----------------------------------------------|------------|
| <b>10 - Field Replaceable Units</b>           | <b>121</b> |
| <b>11 - Disassembly and Repair</b>            | <b>127</b> |
| Tool list                                     | 127        |
| Key Pad Disassembly                           | 132        |
| SpO <sub>2</sub> Circuit Board Disassembly    | 133        |
| Masimo                                        | 133        |
| Nellcor                                       | 134        |
| Blood Pressure Circuit Board Disassembly      | 136        |
| LCD Inverter Ballast Board Disassembly        | 137        |
| LCD Disassembly                               | 137        |
| Pump Disassembly                              | 138        |
| RS232-Communication Circuit Board Disassembly | 138        |
| Fan and Power Circuit Board Disassembly       | 139        |
| Thermometry Circuit Board Disassembly         | 139        |
| SureTemp Plus                                 | 139        |
| Braun ThermoScan PRO 4000                     | 141        |
| <b>12 - Repair Test Specifications</b>        | <b>143</b> |
| General Unit Test                             | 143        |
| A-D Noise Test                                | 143        |
| Leak Test                                     | 143        |
| Inflation Test                                | 143        |
| Dump Test                                     | 143        |
| Pneumatic Accuracy Test                       | 143        |
| Valve Control Test                            | 144        |
| Voltage Calibration                           | 144        |
| Blank Mode Current Test                       | 144        |
| Back Light (Idle) Current Test                | 144        |
| Valve/Pump Mode Current Test                  | 144        |
| Interface Test                                | 144        |
| Temperature Option Requirements               | 144        |
| Accuracy Test                                 | 144        |
| Temperature Probe Test                        | 145        |
| SpO <sub>2</sub> Option Requirements          | 145        |
| SpO <sub>2</sub> Functional Test              | 145        |
| SpO <sub>2</sub> Mode Current Test            | 145        |
| Fail Safe Test                                | 145        |
| Over Pressure Test                            | 145        |
| Over 15 mmHg                                  | 145        |
| Interconnect Diagram 1                        | 147        |
| Interconnect Diagram 2                        | 148        |



# 1 Safety Summary

---

## Introduction

All users of the Spot Vital Signs LXi must read and understand the safety summary, and all additional specific warnings and cautions located throughout the documentation.

## Warnings and Cautions

Familiarize all operating personnel with the general safety information in this summary. Specific warnings and cautions are also found throughout this manual.

### General Warnings

A warning statement in this manual identifies a condition or practice, which if not corrected or discontinued immediately, could lead to patient injury, illness, or death.

These warnings pertain entirely to Spot Vital Signs LXi.



**WARNING** The information in this manual is a comprehensive guide to the operation of Spot Vital Signs LXi. For best results, read this manual thoroughly before using the device.

**WARNING** Spot Vital Signs LXi is designed for medical clinician use. Although this manual may illustrate medical spot-check techniques, only a trained clinician who knows how to take and interpret a patient's vital signs should use this device.

**WARNING** Spot Vital Signs LXi is not intended for use in environments that are without health care practitioner supervision.

**WARNING** Spot Vital Signs LXi is not intended for continuous monitoring. **Do not leave the device unattended while taking measurements on a patient.**

**WARNING** To ensure data integrity, save readings and clear the Spot Vital Signs LXi display between patients.

**WARNING** The Spot Vital Signs LXi is not defibrillator proof.

**WARNING** Spot Vital Signs LXi is not intended for use during patient transport.

**WARNING** This device is not suitable for use in the presence of a flammable anesthetic mixture with air or oxygen or nitrous oxide. An explosion may result.



**WARNING** To ensure patient safety, use only accessories and supplies (i.e., cuffs, hoses, temperature probes, SpO<sub>2</sub> sensors, etc.) recommended for or supplied with Spot Vital Signs LXi. Using unapproved accessories with Spot Vital Signs LXi can affect patient and/or operator safety.

**WARNING** Prevent fluid from entering any connectors on the device. Should this occur, dry the connectors with warm air. Check the accuracy of all operating functions.

**WARNING** Every three months, inspect the blood pressure cuff, SpO<sub>2</sub> cable, and other accessories for fraying or other damage. Replace as necessary.

**WARNING** Do not use Spot Vital Signs LXi on patients who are on heart/lung machines.

**WARNING** Electric shock hazard. There are no user-serviceable parts inside Spot Vital Signs LXi other than battery replacement (see “Battery Replacement” on page 46). Only an individual specifically trained and approved for the repair and/or verification of this specific Welch Allyn product may perform maintenance procedures specifically described in this manual. For service, refer the device to an Authorized Service Center.

**WARNING** This device is not intended for hand-held use during operation.

**WARNING** Do not autoclave.

**WARNING** This device complies with current required standards for electromagnetic interference and should not present problems to other equipment or be affected by other devices. As a precaution, avoid using this device in close proximity to other equipment.

**WARNING** Welch Allyn is not responsible for the integrity of any mounting installation. Welch Allyn recommends that the customer contact their Biomedical Engineering Department or maintenance service to ensure professional installation for safety and reliability of any mounting accessory.

**WARNING** The Spot Vital Signs LXi consists of high-quality precision parts. Protect it from severe impact and shock. A qualified service technician must check any Spot Vital Signs LXi that is dropped or damaged for proper operation prior to further use. Do not use the Spot Vital Signs LXi if you notice any signs of damage. Contact the Welch Allyn Technical Support Department for assistance.

**WARNING** Do not use an SpO<sub>2</sub> finger clip sensor and a blood pressure cuff simultaneously on the same limb. Doing so may result in inaccurate pulse rate and perfusion readings.

**WARNING** All signal input and output (I/O) connectors are intended for connection of only devices complying with IEC 60601-1, or other IEC standards (for example, IEC 60950) as appropriate to the device. Connecting additional devices to the Spot LXi may increase leakage currents. To maintain operator and patient safety, it is necessary to consider the requirements of IEC 60601-1-1.

**WARNING** For proper patient electrical isolation, use only a Welch Allyn power supply (4500-101A) to charge Spot Vital Signs LXi.

**WARNING** Do not use the manufacturer supplied charger for the weight scale while it is attached to Spot LXi.

## Blood Pressure Warnings

These warnings pertain to the Spot Vital Signs LXi blood pressure feature.



**WARNING** Spot Vital Signs LXi is not intended to measure BLOOD PRESSURE on neonatal patients. The AAMI SP10:2002 standard defines neonates as children 28 days or less of age if born at term (37 weeks gestation or more); otherwise up to 44 gestational weeks.

**WARNING** To ensure pediatric blood pressure accuracy and safety, the Child Reusable Two-Piece Blood Pressure Cuff (4500-01), Infant Durable One-Piece Cuff (5082-82-4MQ), and the Infant Disposable One-Piece Cuff (5082-92-4MQ) are the smallest cuffs approved for use with young children and infants. The child's arm must fit within the range markings on the cuff.

**WARNING** Avoid compression of the blood pressure hose or cuff tubing of Spot Vital Signs LXi. This may cause system errors to occur in the device.

**WARNING** Patients who are experiencing moderate to severe arrhythmias may give inaccurate blood pressure measurements.

**WARNING** Spot Vital Signs LXi does not operate effectively on patients who are experiencing convulsions or tremors.

**WARNING** Use only Welch Allyn blood pressure cuffs and/or hoses. Using other manufacturers' blood pressure cuffs and/or hoses may produce inaccurate blood pressure readings.

**WARNING** When several blood pressure measurements are taken on the same patient, regularly check the cuff site and extremity for possible ischemia, purpura, and/or neuropathy.

**WARNING** Do not place the cuff on any extremity that is used for intravenous infusions or any area where circulation is compromised.

**WARNING** Excessive cuff tightness may cause venous congestion and discoloration of the limb.

**WARNING** Wrapping the cuff too loosely (preventing proper inflation) may result in errors.

**WARNING** Do not change the connector(s) on the blood pressure cuff tubing of this device to luer type. Luer type connectors are commonly used in intravenous infusion systems. Using the luer connectors on blood pressure cuff tubing creates the risk that the blood pressure tubing could be mistakenly connected to a patient's intravenous line, resulting in the introduction of air into the patient's circulatory system.

## Temperature Warnings

These warnings pertain to the Spot Vital Signs LXi temperature feature.

### SureTemp® Plus

These warnings are specific to the SureTemp Plus thermometer option.



**WARNING** Use only Welch Allyn probe covers. Using other manufacturers' probe covers or no probe cover may produce temperature measurement errors and/or inaccuracy.

**WARNING** Always use a probe cover whenever coming into contact with a patient.

**WARNING** Long-term continuous monitoring beyond three to five minutes is not recommended in any mode.

**WARNING** Oral/axillary probes (blue ejection button at top of probe) and blue removable probe wells are used for taking oral and axillary temperatures only. Rectal probes (red ejection button) and red removable probe wells are used for taking rectal temperatures only. Use of the probe at the wrong site will result in temperature errors. Use of the incorrect removable probe well could result in patient cross-contamination.

**WARNING** The thermometer connectors and probe are not waterproof. Do not immerse or drip fluids on these items. Should this occur, dry the connectors and probe with warm air. Check all functions for proper operation and accuracy.

**WARNING** Do not take an axillary temperature through patient's clothing. Direct probe cover to skin contact is required.

**WARNING** Do not autoclave.

**WARNING** Use Welch Allyn single-use disposable probe covers to limit patient cross-contamination.

**WARNING** Incorrect insertion of probe can cause bowel perforation.

**WARNING** Washing hands greatly reduces the risk of cross-contamination and nosocomial infection.

**WARNING** To ensure optimal accuracy, always confirm that the correct mode is selected.

## Braun ThermoScan PRO 4000

These warnings are specific to the Braun ThermoScan PRO 4000 thermometer option.



**WARNING** Keep the probe window clean, dry, and undamaged at all times to ensure accurate measurements. To protect the probe window, always keep the thermometer in the storage cover while transporting or when not in use.

**WARNING** Only use Braun ThermoScan probe covers with this thermometer. Using other manufacturer's probe covers or no probe cover may produce temperature measurement errors and/or inaccuracies. If the thermometer is used without a probe cover attached, clean the lens (see "Braun ThermoScan PRO 4000 Thermometer" on page 44).

**WARNING** Do not autoclave.

**WARNING** The thermometer is not waterproof. Do not immerse or drip fluids on it. Should this occur, dry the thermometer with warm air. Check for proper operation and accuracy.

## SpO<sub>2</sub> Warnings

These warnings pertain to the Spot Vital Signs LXi SpO<sub>2</sub> feature.



**WARNING** Only use Spot Vital Signs LXi with Masimo or Nellcor SpO<sub>2</sub> option with Masimo or Nellcor brand sensors and accessories, respectively. Using the wrong or unapproved sensors or cables may cause improper performance.

**WARNING** The SpO<sub>2</sub> sensor and extension cables are intended for use only for pulse oximetry measurements. Do not attempt to connect these cables to a PC or any similar device.

**WARNING** Before using, carefully read the sensor Directions for Use, including all warnings, cautions, and instructions.

**WARNING** Do not use a damaged sensor or pulse oximetry cable or a sensor with exposed optical components.

**WARNING** Incorrect application or a long duration of use of an SpO<sub>2</sub> sensor may cause tissue damage. Inspect the sensor site periodically as directed in the sensor's Directions for Use.

**WARNING** Certain ambient environmental conditions, sensor application errors, and certain patient conditions may affect SpO<sub>2</sub> readings and pulse signal.

**WARNING** Do not immerse the sensor or patient cables in water, solvents, or cleaning solutions (the sensors and connections are not waterproof). Do not use irradiation, steam, or ethylene oxide for sterilization.

**WARNING** The SpO<sub>2</sub> in the Spot Vital Signs LXi device is not intended for use as an apnea monitor.

**WARNING** Consider the SpO<sub>2</sub> an early warning device. As a trend toward patient deoxygenation is indicated, use laboratory instruments to analyze blood samples to completely understand the patient's condition.



**WARNING** Tissue damage can be caused by incorrect application or duration of use of a Nellcor OxiMax sensor. Inspect the sensor site as directed in the sensor Directions for Use.

**WARNING** Do not use the sensors during magnetic resonance imaging (MRI) scanning. Induced current could potentially cause burns. The MS board pulse oximeter may affect the MRI image, and the MRI unit may affect the accuracy of the oximetry measurements.

**WARNING** Carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.

**WARNING** Failure to cover the Nellcor OxiMax sensor site with opaque material in high ambient light conditions may result in inaccurate measurements.

## General Cautions

A caution statement in this manual identifies a condition or practice, which if not corrected or discontinued immediately, could lead to equipment failure, equipment damage, or data loss.

These cautions pertain to the entire Spot Vital Signs LXi device.



**Caution** If the accuracy of any measurement is in question, check the patient's vital sign(s) with an alternate method and then check to verify the device is functioning properly.

**Caution** Place the device on a secure surface or use one of the optional mounting accessories.

**Caution** Do not place fluids on or near the device.

**Caution** It is recommended that the device is used within stated operating temperature ranges (see "Environmental" on page 93). The device will not meet its performance specifications if used outside these temperatures ranges.

**Caution** Always unplug the AC power transformer from the outlet before moving the mobile stand to a new location.

**Caution** The basket has a three-pound weight limit. Take care not to exceed this limit.

**Caution** When connecting a Health o meter<sup>®</sup> weight scale, only operate the scale with battery power. Remove and discard the scale's AC power cord.

**Caution** When using the Health o meter scale, remove the two hex nuts on the RS-232 cable, as supplied, before screwing the cable into the scale base.

## Blood Pressure Cautions

These cautions pertain to the Spot Vital Signs LXi blood pressure feature.



**Caution** Minimize extremity and cuff motion during blood pressure readings.

**Caution** If the blood pressure cuff is not at heart level, note the difference in reading due to the hydrostatic effect. Add the value of 1.80 mmHg (.2 kPa) to the displayed reading for every inch (2.5 cm) above heart level. Subtract the value of 1.80 mmHg (.2 kPa) from the displayed reading for every inch (2.5 cm) below heart level.

**Caution** Proper blood pressure cuff size and placement is essential to the accuracy of the blood pressure determination. See Reusable Two-Piece Cuff Measurements or Durable One-Piece Cuff Measurements of the Directions for Use for sizing information.

**Caution** The position and physiologic condition of the subject can affect a blood pressure reading.

## Temperature Cautions

These cautions pertain to the Spot Vital Signs LXi temperature feature.



**Caution** The SureTemp Plus feature only operates with the probe well in place.

**Caution** Do not use alkaline batteries in the Braun ThermoScan PRO 4000.

**Caution** Biting the probe tip may result in damage to the probe.

## SpO<sub>2</sub> Cautions

These cautions pertain to the Spot Vital Signs LXi SpO<sub>2</sub> feature.



**Caution** The pulse oximeter is calibrated to determine the percentage of arterial oxygen saturation of functional hemoglobin. Significant levels of dysfunctional hemoglobin such as carboxyhemoglobin or methemoglobin may affect the accuracy of the measurement.

**Caution** Some sensors may not be appropriate for a particular patient. If at least 10 seconds of perfusion pulses cannot be observed for a given sensor, change sensor location or sensor type for perfusion to resume.

**Caution** Physiological conditions, medical procedures, or external agents that may interfere with the pulse oximeter's ability to detect and display measurements include dysfunctional hemoglobin, arterial dyes, low perfusion, dark pigment, and externally applied coloring agents such as nail polish, dye, or pigmented cream.

**Caution** When selecting a sensor, consider the patient's weight and activity level, the adequacy of perfusion, the available sensor sites, the need for sterility, and the anticipated duration of monitoring.

## Electrostatic Discharge (ESD)



Electrostatic discharge is a sudden current flowing from a charged object to another object or to ground. Electrostatic charges can accumulate on common items such as foam drinking cups, cellophane tape, synthetic clothing, untreated foam packaging material, and untreated plastic bags and work folders, to name only a few.

Electronic components and assemblies, if not properly protected against ESD, can be permanently damaged or destroyed when near or in contact with electrostatically charged objects. When you handle components or assemblies that are not in protective bags and you are not sure whether they are static-sensitive, assume that they are static-sensitive and handle them accordingly.

- Perform all service procedures in a static-protected environment. Always use techniques and equipment designed to protect personnel and equipment from electrostatic discharge.
- Remove static-sensitive components and assemblies from their static-shielding bags only at static-safe workstations - a properly grounded table and grounded floor mat - and only when you are wearing a grounded wrist strap (with a resistor of at least 1 megohm in series) or other grounding device.
- Use only grounded tools when inserting, adjusting, or removing static-sensitive components and assemblies.
- Remove or insert static-sensitive components and assemblies only with monitor power turned off.
- Insert and seal static-sensitive components and assemblies into their original static-shielding bags before removing them from static-protected areas.

Always test your ground strap, bench mat, conductive work surface, and ground cord before removing components and assemblies from their protective bags and before beginning any disassembly or assembly procedures.



# Symbols

The following symbols are associated with the Spot Vital Signs LXi.

## Safety Symbols



Identifies information within the manual to avoid injury.



Identifies information within the manual to avoid equipment failure.



Caution: consult accompanying documents



Internally Powered, Lead Acid Battery



Handle with Care



Transport Temperature



Storage Humidity



Recycle



Class II Equipment



Equipment is not protected against the ingress of liquid.



Type BF Equipment



On/Off



Recycle the product separate from other disposables, see "Product Disposal" on page 48.



Non-ionizing radiation (RF transmitter)

Mode of Operation: Continuous



DC Power In

## Button Symbols



Navigation Buttons



Select



Blood Pressure



Memory



Power On/Off

## Connection Symbols



USB Connection



Serial Port Connection


## Agency Symbols

CONFORMS TO:  
UL STD 60601-1

IEC 60601-1



The CE mark on this product indicates that it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.



Regulatory Affairs Representative  
Welch Allyn Limited  
Navan Business Park  
Dublin Road  
Navan, County Meath, Republic of Ireland

# 2 Overview

---

## Purpose and Scope

The Spot Vital Signs LXi Service Manual is intended as a reference for maintenance and repair to the field replaceable unit (FRU) level and are listed on [page 121](#).

**Note** Opening the Spot Vital Signs device case requires recalibration of specific functions. This includes any internal components that were replaced. Calibration requires the use of the Welch Allyn Spot LXi Tool.

This manual provides the technical qualified service person with troubleshooting information, repair procedures, and calibration and performance verification instructions. A technical overview of the Spot Vital Signs LXi subsystems is provided as an introduction to the device's circuitry and pneumatics.

This manual is intended for the technical qualified service personnel. Service training classes on Welch Allyn's products are available. Contact Welch Allyn Technical Support for information.

## Other Applicable Documents

The Spot Vital Signs LXi Directions for Use manual is also available. Refer to this document for information other than maintenance and repair.

Welch Allyn 9600 Plus Calibration Tester Directions for Use - for all models.


Braun ThermoScan PRO 4000 User's Guide - for models 450E0, 45NE0, 45ME0.

Masimo Directions for Use - for models 45MT0, 45ME0.

Nellcor Directions for Use - for models 45NT0, 45NE0.


# Miscellaneous Mounting Accessories

## Wall Mount Kit




| PARTS LIST |                  |     |           |
|------------|------------------|-----|-----------|
| ITEM       | DESCRIPTION      | QTY | P/N       |
| 1          | BASKET           | 1   | MIS-1574M |
| 2          | AC ADAPTER STRAP | 1   | STL-1534  |
| 3          | ADAPTER BLOCK    | 1   | PLS-1497  |
| 4          | BLOCK CAP        | 1   | STL-1534  |
| 5          | SCREW            | 2   | SCR-1464  |

UTILITY BASKET  
12.00" X 8.12" X 7.00"



WALL MOUNT BASKET WITH GRAY EPOXY FINISH. SPACE SAVING UNDER MOUNT AC ADAPTER HOLDER.  
BULK SHIPPING PACK OF 8 UNITS. SHIPPING BOX IS 27.00" X 18.00" X 18.00" SHIPPING WEIGHT IS 15 LBS.  
ITEM 1 TO BE PLACED IN A POLY BAG, ITEMS 2-5 TO BE PLACED IN A POLY BAG.


## Mobile Stand Kit



5-LEG, 22" DIAMETER LOW CENTER-OF-GRAVITY DETACHABLE BASE, WITH 3" SOFT RUBBER CASTERS, 2 EA. LOCKING. GRAY EPOXY FINISH. DETACHABLE POLE ASSEMBLY WITH GRAY BASKET AND HANDLE, WITH TILT MOUNT, AND WITH DEVICE MOUNTING PLATE. SHIPPING BOX IS 40.00" X 21.50" X 12.50". SHIPPING WEIGHT IS 31 LBS.


NOTE: ALL DIMENSIONS ARE IN INCHES.

| PARTS LIST |                      |      |           |
|------------|----------------------|------|-----------|
| ITEM       | DESCRIPTION          | QTY  | P/N       |
| 1          | STEEL BASE           | 1    | STL-1005  |
| 2          | STUD                 | 5    | SCR-1258  |
| 3          | CASTER               | 3    | CAS-1151  |
| 4          | ESNA NUT             | 5    | SCR-1257  |
| 5          | WASHER               | 5    | SCR-1341  |
| 6          | LOCKING CASTER       | 2    | CAS-1163  |
| 7          | WASHER               | 1    | SCR-1274  |
| 8          | HEX BOLT             | 1    | SCR-1349  |
| 9          | TAPERED PIN          | 1    | STL-1032  |
| 10         | BASE POLE            | 1    | TUB-1268  |
| 11         | BUSHING              | 1    | PLS-1436  |
| 12         | BASKET               | 1    | MIS-1574M |
| 13         | MOUNT                | 1SET | PLS-1496  |
| 14         | TUBE                 | 1    | TUB-1202  |
| 15         | END CAPS             | 2    | MIS-1461  |
| 16         | SCREW                | 1    | SCR-1567  |
| 17         | SWIVEL MOUNT         | 1    | PLS-1470  |
| 18         | SCREWS               | 3    | SCR-1356  |
| 19         | MOUNT PLATE          | 1    | PLS-1502  |
| 20         | KNOB                 | 1    | KNB-1339  |
| 21         | RETAINING WASHER     | 1    | SCR-1575  |
| 22         | SCREW                | 2    | SCR-1295  |
| 23         | WASHER               | 2    | SCR-1402  |
| 24         | ACORN NUT            | 2    | SCR-1372  |
| 25         | SCREW                | 2    | SCR-1559  |
| 26         | WASHER               | 2    | SCR-1402  |
| 27         | ESNA NUT             | 2    | SCR-1371  |
| 28         | AC ADAPTER STRAP     | 1    | STL-1534  |
| 29         | ADAPTER BLOCK        | 1    | PLS-1497  |
| 30         | BLOCK CAP            | 1    | STL-1534  |
| 31         | SCREW                | 2    | SCR-1464  |
| 32         | MFG LABEL            | 1    | LBL-1526  |
| NOT SHOWN  | TAMPER PROOF HEX KEY | 1    | SCR-1320  |
|            | TAMPER PROOF SCREW   | 1    | SCR-1574  |
|            | CARTON               | 1    | PAC-1530  |
|            | INSERT               | 1    | PAC-1535  |
|            | BASE BOX             | 1    | PAC-1536  |
|            | BASKET BOX           | 1    | PAC-1537  |
| PACKAGING  | STYROFOAM INSERTS    | 2    | PAC-1538  |
|            | INSTRUCTION SHEET    | 1    | LBL-1948  |




TORQUE ITEM 27 TO 8 - 10 in lbs.  
BAG ITEMS 7,8 AND 28 - 31, ALONG WITH TAMPER PROOF KEY AND SCREW.

## Radio Mounting Accessory





## Barcode Reader Mounting Accessory



# Warranty

## Spot Vital Signs LXi

Welch Allyn warrants that the Spot Vital Signs LXi Radio products meet the labeled specifications of the products and will be free from defects in materials and workmanship that occur within one year after the date of purchase.

The date of purchase is: 1) the date specified in our records if you purchased the Product directly from us, 2) the date specified in the warranty registration card that we ask you to send to us, or 3) the date of purchase of product from the authorized Welch Allyn distributor as documented from a receipt from said distributor.

This warranty does not cover damage caused by: 1) handling during shipping, 2) use or maintenance contrary to labeled instructions, 3) alteration or repair by anyone not authorized by Welch Allyn, and 4) accidents.

If a Product or accessory covered by this warranty is determined to be defective because of defective materials, components, or workmanship, and the warranty claim is made within the warranty period described above, Welch Allyn will, at its discretion, repair or replace the defective Product or accessory free of charge.

You must obtain a service notification number from Welch Allyn to return your Product before you send it to Welch Allyn's designated service center for repair. Contact Welch Allyn Technical Support.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. WELCH ALLYN'S OBLIGATION UNDER THIS WARRANTY IS LIMITED TO REPAIR OR REPLACEMENT OF PRODUCTS CONTAINING A DEFECT. WELCH ALLYN IS NOT RESPONSIBLE FOR ANY INDIRECT OR CONSEQUENTIAL DAMAGES RESULTING FROM A PRODUCT DEFECT COVERED BY THE WARRANTY.

## Accessories

Refer to the manufacturer's Directions for Use for the Masimo and Nellcor finger sensor and cable warranty.

Refer to the Directions for Use provided with the Welch Allyn Blood Pressure Cuff for warranty information.

The SureTemp Plus probe is covered by a one-year warranty and the SureTemp Plus probe well is covered by a 90-day warranty against original defects in material and workmanship. Probe covers are intended for single-use only.

The Braun ThermoScan PRO 4000 is covered by a three-year warranty against original defects in material or workmanship.

The printer is covered by a one-year warranty against original defects in material or workmanship.

The barcode scanner is covered by a five-year warranty against original defects in material or workmanship.

# Service



**Caution** Unauthorized repairs will void the warranty.

All repairs on products under warranty must be completed by Welch Allyn or Welch Allyn certified service personnel in accordance with warranty guidelines. See the Welch Allyn Partners In Care™ web site (<http://www.welchallyn.com/promotions/services/default.htm>) for a description of the different levels of services offered.

## User Serviceable Parts

For user serviceable parts available see “Supplies and Accessories” on page 113.

## Technical Assistance

If you have an equipment problem that you cannot resolve, call the Welch Allyn Service Center nearest you for assistance. Technical service telephone support is available on normal business days.

If you are advised to return a product to Welch Allyn for repair or routine maintenance, schedule the repair with the service center nearest you.

Before returning a product for repair, you must obtain authorization from Welch Allyn. Service personnel will give you a Service Notification number. Please note this number on the outside of your shipping box. Returns without a Service Notification number will not be accepted for delivery.

## Field Replacement Units

See “Field Replaceable Units” on page 121.

## Service Loaners

Service loaners or exchanges may be provided depending on the specific Partners In Care™ agreement in place for the device. Loaners for products repaired while under the original warranty, may be provided free of charge and are shipped within 48 hours of notification of need.

For service repairs outside of warranty or outside of service agreement are available for a nominal charge and shipment is subject to availability. Loaners are shipped pre-paid; however, this charge is added to the service charges.

## Service Intervals

Calibrate the Spot Vital Signs LXi on an annual basis for blood pressure, temperature, and SpO2 accuracy.


# Spot Vital Signs LXi Configurations


**Table 1. Available Versions of Spot Vital Signs LXi**

| REF   | Description                                                                    |
|-------|--------------------------------------------------------------------------------|
| 450T0 | SureBP with SureTemp Plus Thermometer                                          |
| 450E0 | SureBP with Braun ThermoScan PRO 4000 Thermometer                              |
| 45MT0 | SureBP with Masimo SpO <sub>2</sub> and SureTemp Plus Thermometer              |
| 45ME0 | SureBP with Masimo SpO <sub>2</sub> and Braun ThermoScan PRO 4000 Thermometer  |
| 45NT0 | SureBP with Nellcor SpO <sub>2</sub> and SureTemp Plus Thermometer             |
| 45NE0 | SureBP with Nellcor SpO <sub>2</sub> and Braun ThermoScan PRO 4000 Thermometer |

## Controls

**Figure 1. Spot Vital Signs LXi Front Panel with SureTemp Plus Thermometer**




**Figure 2. Spot Vital Signs LXi Front Panel with Braun ThermoScan PRO 4000 Thermometer**

# Display Window

**Note** Before using Spot Vital Signs LXi for the first time, you must program it using the initial configuration screen.

1. Press the **Power** button. The display window shows the initial configuration screen.

**Figure 3. Initial Configuration Screen**



The image shows a rectangular display window with a black border. Inside, the text is as follows: At the top, 'Change Local Defaults Menu' is centered. Below it, 'Model No.: 4 NTO' and 'Serial No. 200 040004' are centered. Then, a list of settings is shown on the left: 'Language:', 'BP Units:', 'Temperature Units:', 'eight Units:', 'Weight Units:', 'Date Format:', 'Date:', 'Time Format:', and 'Time:'. At the bottom right, there is a horizontal line with the word 'Select' centered below it.

Change Local Defaults Menu

Model No.: 4 NTO

Serial No. 200 040004

Language:

BP Units:

Temperature Units:

eight Units:

Weight Units:

Date Format:

Date:

Time Format:

Time:


\_\_\_\_\_

Select

2. Use the **Select** button to access the options and accept the entries, and use the **Navigation** buttons to move through the menu.
3. The word "Exit" appears at the bottom of the list after you have programmed all items in the menu. You must program all items before you can start to use the device.
4. Scroll to Exit and press the **Select** button to save the entries.

The liquid crystal display (LCD) may indicate any of the following: systolic blood pressure (mmHg or kPa), diastolic blood pressure (mmHg or kPa), MAP (mmHg or kPa), temperature (°F or °C), temperature mode, pulse rate, pulse signal level, SpO<sub>2</sub> percent, department location, date, time, record number, height (in or cm), weight (lb or kg), respiration rate, pain level, connectivity signal strength, and battery charge level.


**Figure 4. Display Window**



# Connections

Use the following instructions to connect the blood pressure hose, thermometer probe, and optional attachments to the Spot Vital Signs LXi.

**Figure 5. Spot Vital Signs LXi Side and Rear Panel Connections**






## Braun ThermoScan PRO 4000 Lock

Press the lock tab toward the Braun ThermoScan PRO 4000 thermometer until it clicks. To release the housing, insert the lock release pin into the lock release hole until the lock tab snaps back.

**Figure 6. Spot Vital Signs LXi with Braun ThermoScan PRO 4000 Lock**



## Blood Pressure Hose and Cuff

Identify and have available the Spot Vital Signs LXi, blood pressure cuff, and the blood pressure hose.

1. Inspect the blood pressure hose; notice that one end has a single, gray connector fitting and the other end has two white fittings.
2. Squeeze the side tabs on the gray connector and completely push the blood pressure hose connector into the blood pressure hose connection port until it clicks into place (see [Figure 5](#) on page 24).
3. Twist the white connectors on the blood pressure hose and cuff connectors together.

## Thermometer

Spot Vital Signs LXi is available with either the SureTemp Plus thermometer or the Braun ThermoScan PRO 4000 thermometer.

### SureTemp Plus

SureTemp Plus is available with two probes and matching wells; one for oral/axillary temperatures (blue ejection button and probe well) and one for rectal temperatures (red ejection button and probe well).



**WARNING** Always use a probe cover whenever coming into contact with a patient.



**Caution** The SureTemp Plus feature only operates with the probe well in place.

1. Align the probe well with the tabs facing up and down into the round opening of the SureTemp Plus housing on the right side of Spot Vital Signs LXi. Push it into place.
2. Align the temperature probe connector with the SureTemp Plus thermometer connection port on the back of the Spot Vital Signs LXi (see [Figure 5](#) on page 24). You can only insert the connector into the port one way.
3. Press the tab on the connector and push it until it clicks into place.
4. Insert the temperature probe into the probe well.

### Braun ThermoScan PRO 4000



**Caution** Do not use alkaline batteries in the Braun ThermoScan PRO 4000 thermometer. Use the rechargeable battery pack supplied by Welch Allyn.

1. Open the package of rechargeable batteries and follow the installation guide provided.
2. Open the box of probe covers as directed on the box and slide the box into the metal guides toward the back of the thermometer housing with the opening at the top and perforation facing forward.
3. Hold the Braun thermometer at a 45° angle then insert the probe and the top of the thermometer into the housing.
4. Lower the bottom portion of the thermometer into the housing until it snaps into place. If you do not properly seat the thermometer, it could fall out of the holder and become damaged.
5. Slide the thermometer housing into the thermometer slot on the right side of the Spot Vital Signs LXi device.
6. Push the lock tab forward to prevent the thermometer housing from falling out of Spot Vital Signs LXi (see [Figure 6](#) on page 25).

To release the lock, insert the lock release pin into the lock release hole.

## SpO<sub>2</sub> Sensor


Spot Vital Signs LXi is available with a wide variety of SpO<sub>2</sub> sensors and ships with a reusable finger clip sensor.

1. Align the shape and pin configuration of the extension cable connector to the SpO<sub>2</sub> cable connection port on the left side of the Spot Vital Signs LXi device.
2. Push the connector firmly into the SpO<sub>2</sub> cable connection port until you hear it click into place (see [Figure 5](#) on page 24).
3. Align the opposite end of the extension cable to the sensor cable connector and firmly push them together.



**WARNING** Use only Masimo or Nellcor SpO<sub>2</sub> sensors and accessories with the Spot Vital Signs LXi with Masimo or Nellcor configurations, respectively. Using the wrong or unapproved sensors or cables may cause improper performance.

## Radio



The Spot Vital Signs LXi identifies which serial port to use after activating the Radio in the internal configuration mode.

Welch Allyn offers two different radios:

- 802.11b radio and
- 802.11a/b/g radio

Only use Welch Allyn Radios

## Quick Reference Card

Attach the Quick Reference Card to the Spot Vital Signs LXi handle, mobile stand, or wall mount using the supplied plastic cable tie.

## AC Power Transformer

**Note** To assure proper electrical isolation use only the Welch Allyn SPOT Vital Signs AC power transformer/charger.

The operator can use the Spot Vital Signs LXi with AC or battery power (after charging the battery).

1. Insert the round transformer connector into the AC power connection port on the back of the Spot Vital Signs LXi (see [Figure 5](#) on page 24).
2. Insert the line cord into the line connector on the transformer then plug the power cord on the transformer into the AC main power source to charge the battery.


## Printer

Once the printer is enabled in the Internal Configuration Mode, you can use plain or label paper. If you use label paper, you must clean the printer heads after every 10 rolls. A head cleaner kit is provided with each package of label paper available from Welch Allyn.

### To connect the printer to the Spot LXi:

1. Connect the the large end of the cable (405210) to the Spot LXi.
2. Connect the small end of the cable (405210) to the Spot LXi printer.

Connect to Spot  
LXi



Connect to Spot  
LXi printer

### To charge the Spot LXi printer:

1. Remove the guard on the power supply body (713703 or BC1300M) and attach the appropriate power supply plug/adaptor.

2. Connect the Spot LXi printer power supply to the Spot LXi printer.
3. Charge the Spot LXi printer for at least 9 hours before initial use.

**To enable the Spot LXi printer:**

1. Connect the Spot LXi printer to the Spot LXi.
2. Enable the printer in the Internal Configuration Mode of the Spot LXi (see [“Internal Configuration Mode” on page 32](#)).

## Battery

Charge the Spot LXi battery and printer battery (if using) for 9 hours before initial use. Charge the device an additional hour if it includes a Braun ThermoScan PRO 4000 thermometer.

**Note** While charging the Spot LXi, it is acceptable to charge the Spot LXi printer with the power supply (713703 or BC1300M).

While Spot Vital Signs LXi is charging, the charging LED (~) flashes and the battery level indicator segments on the display continuously sequence. When the battery is fully charged, the charging LED stops flashing and the battery level indicator will stop sequencing.

If the device includes a Braun ThermoScan PRO 4000 thermometer, the charging LED below the thermometer will illuminate orange as it is charging. When the battery is fully charged, the LED will power off.

If the Spot LXi has a mounted, external printer that is connected via the cable, the charging LED on the printer flashes orange and green when the Spot LXi is plugged into the AC main power source. If the Spot LXi has a mounted, external printer that is not connected via the cable, the charging LED on the printer will not flash.

**Note** There is no hazard associated with leaving the battery in the device, even if the device is not used for long periods of time.



# 3

## Functional Overview

---

This functional verification procedure helps to confirm the proper operation of the Spot Vital Signs LXi and options. This procedure supports the requirements of routine preventative maintenance. It is not necessary to disassemble the Spot Vital Signs LXi to perform this procedure.

For the calibration procedures, see [“Calibration and Verification”](#) on page 49. If the Spot Vital Signs LXi fails certain functional tests or a circuit board is replaced, the device may require calibration. Always perform this functional verification procedure after performing any calibration.

### Self Test

The Spot Vital Signs LXi performs a self-test each time the device is powered on. Press the **Power** button to turn the device on or off. Upon each power up, the display turns on, two audible beeps occur, and the Spot Vital Signs LXi displays the model and serial numbers. If the internal self-check is successful, the display shows its normal functions (see [Figure 4](#) on page 23) with all values blank, and the device is ready for operation. If the self-check fails, an error code is shown in the Navigation Window (see “Error Codes” on page 99).

Spot Vital Signs LXi automatically powers off when not used for 30 minutes.


# Internal Configuration Mode

You can change several device operating parameters in the Internal Configuration Mode. When changed, these settings become the default power-up settings. You will also see non-changeable device configurations for technical service purposes.

## To Enter the Internal Configuration Mode:

1. Turn the Spot Vital Signs LXi off.
2. Press and hold the **Select** and **Power** buttons for 5 seconds. The device enters the Internal Configuration Mode and the Configuration Menu screen appears on the display.

**Figure 7. Internal Configuration Mode Menu**



3. Use the **Navigation** buttons to move through the menu options and then press the **Select** button to access the options or accept a change. See the following tables for descriptions of the menu options.
4. Press the **Power** button to exit the Internal Configuration Mode.



**Table 2. Configuration Menu Options**

| Setting             | Description                                                                                                                                                                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Version Numbers*    | Displays the software and hardware version numbers in the Spot Vital Signs LXi device.                                                                                                                                                                           |
| Battery*            | Displays the battery level.                                                                                                                                                                                                                                      |
| Location Identifier | Allows the entry of the device's location (e.g., the department name). Follow the display prompts and use the <b>Navigation</b> and <b>Select</b> buttons to enter up to 10 characters.                                                                          |
| Date/Time           | Changes the date and time formats or updates the actual date and time. See <a href="#">Table 3, "Date/Time Menu Options"</a> for available settings.                                                                                                             |
| Blood Pressure      | Changes the blood pressure options. See <a href="#">Table 4, "Blood Pressure Menu Options"</a> for available settings.                                                                                                                                           |
| Temperature         | Changes the temperature options. See <a href="#">Table 5, "Temperature Menu Options"</a> for available settings.                                                                                                                                                 |
| Contrast            | Changes the Display Contrast options. Use the left/right navigation buttons to adjust.                                                                                                                                                                           |
| Manual Parameters   | Changes the manual parameters defaults. See <a href="#">Table 6, "Manual Parameters Menu Options"</a> for available settings.                                                                                                                                    |
| External Devices    | Enables or disables available external devices. See <a href="#">Table 7, "External Devices Menu Options"</a> for available settings.                                                                                                                             |
| Save Readings       | Saves the current patient reading at a preselected time interval or upon request. See <a href="#">Table 8, "Save Readings Menu Options"</a> for available settings.                                                                                              |
| Buttons Lockout     | Secures Spot Vital Signs LXi so unauthorized people cannot use the device or access data without enacting the proper key sequence. To override the buttons lockout feature, simultaneously press the <b>Left Navigation</b> button and the <b>Select</b> button. |
| Defaults            | Allows the user to select the default settings for the device and reset the unit to the default settings. See <a href="#">Table 9, "Change Local Defaults Options"</a> for available settings.                                                                   |
| Event Log*          | Displays the recent button presses, errors, measurements, measurement sites, battery state changes, and patient reading send events.                                                                                                                             |

\* Displayed information only; operator cannot change.

**Table 3. Date/Time Menu Options**

| Setting     | Description                                                                                                                                                                                                                                                          |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Date Format | Displays the date in one of the following styles: <ul style="list-style-type: none"> <li>• mm/dd/yyyy example: July 16, 2005 = 07/16/2005</li> <li>• dd/mm/yyyy example: 16 July 2005 = 16/07/2005</li> </ul>                                                        |
| Date        | Changes the date on the Display Window and in patient readings. If a location is entered (see " <a href="#">Location Identifier</a> " in <a href="#">Table 2, "Configuration Menu Options"</a> ), the date will not appear on the Display Window; the location will. |
| Time Format | Displays the time in one of the following styles: <ul style="list-style-type: none"> <li>• 12-hour example: 5:00 PM</li> <li>• 24-hour example: 17:00</li> </ul>                                                                                                     |
| Time        | Changes the time on the Display Window.                                                                                                                                                                                                                              |

**Table 4. Blood Pressure Menu Options**

| Setting                      | Description                                                                                                                                                                                                   |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| BP Calibration Check         | Prepares the Spot Vital Signs LXi for calibration. Only qualified personnel should verify the Spot Vital Signs LXi blood pressure calibration. For more details, see "Blood pressure calibration" on page 62. |
| Blood Pressure Units         | mmHg or kPa.                                                                                                                                                                                                  |
| Mean Arterial Pressure (MAP) | On or off.                                                                                                                                                                                                    |

**Table 5. Temperature Menu Options**

| Setting           | Description                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Temperature Units | Fahrenheit (°F) or Celsius (°C).                                                                                                                                                                                                                                                                          |
| Temperature Mode  | SureTemp Plus models only: Oral, Pediatric Axillary, Adult Axillary, and Last Mode. In Last Mode the device takes the next temperature in the mode in which the previous temperature was measured. Rectal Mode is available only when the rectal probe (red ejection button) and probe well are attached. |

**Table 6. Manual Parameters Menu Options**

| Setting        | Description                                                                                                                               |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Height         | On or off.                                                                                                                                |
| Height Units   | Inches (in) or centimeters (cm).                                                                                                          |
| Height Default | Changes the default patient height displayed in the Navigation Window.                                                                    |
| Weight         | On or off. Even if weight is enabled here, if weight scale is enabled in the External Devices Menu, you cannot manually enter the weight. |
| Weight Units   | Pounds (lb) or kilograms (kg).                                                                                                            |
| Weight Default | Changes the default patient weight displayed in the Navigation Window.                                                                    |
| Respiration    | On or off.                                                                                                                                |
| Pain Level     | On or off.                                                                                                                                |

**Table 7. External Devices Menu Options**

| Setting               | Description                                                                                                                                                                  |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Information System    | On or off. You must enable this option to send patient readings wired or wirelessly.                                                                                         |
| Barcode Patient ID    | On or off. You must enable this option to send patient readings wirelessly.                                                                                                  |
| Barcode Clinician ID  | On or off. You must enable this option to require mandatory sign in of a clinician.                                                                                          |
| Mandatory Sign In     | Yes or no. You must enable this option to require the sign-in of the clinician.                                                                                              |
| Clear on Send/Save:   | Yes or no. You must enable this option to clear the Clinician ID after sending or saving the measurements.                                                                   |
| Clear on Power Cycle: | Yes or No. You must enable this option to clear the Clinician ID after powering the Spot Vital Signs LXi off.                                                                |
| Weight Scale          | On or off. Spot Vital Signs LXi connects to a scale and the weight shows on the display window (see "Spot Vital Signs LXi Technical Overview" on page 37 for scale details). |
| Wireless Module       | None or DPAC. You must enable 802.11 B to send patient readings wirelessly. The radio is available as an accessory only with Welch Allyn Connex.                             |
| Printer               | On or off.                                                                                                                                                                   |
| Printer Paper         | Plain or labels. Only available if the Printer is enabled.                                                                                                                   |

**Table 8. Save Readings Menu Options**

| Setting             | Description                                                                                                                                                                                                                                                                                                                         |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Save Mode           | Manual or automatic. If automatic, Spot Vital Signs LXi saves readings at a preselected time interval. For either option, Spot Vital Signs LXi automatically saves the measured parameters into memory before automatically powering off when not used for 30 minutes.                                                              |
| Auto Save Interval  | Changes the amount of time before automatically saving the current patient reading. Only available if Automatic Save Mode is enabled.                                                                                                                                                                                               |
| Reading Full Action | Auto Overwrite, Ask Overwrite, Do not Overwrite. Spot Vital Signs LXi can save 50 patient readings in memory. Upon reaching reading 51, the device may automatically overwrite reading 1, ask the user if he/she wants to overwrite reading 1, or disable the ability to take another reading until at least one reading is erased. |

**Table 9. Change Local Defaults Options**

| Setting           | Description                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Language          | English, Dansk, Nederlands, Suomi, Français, Deutsch, Italiano, Norsk, Español, Português, Svenska, or Chinese.                                                                                               |
| BP Units          | mmHg or kPa.                                                                                                                                                                                                  |
| Temperature Units | Fahrenheit (°F) or Celsius (°C).                                                                                                                                                                              |
| Height Units      | Inches (in) or centimeters (cm).                                                                                                                                                                              |
| Weight Units      | Pounds (lb) or kilograms (kg).                                                                                                                                                                                |
| Date Format       | Displays the date in one of the following styles: <ul style="list-style-type: none"> <li>• mm/dd/yyyy example: July 16, 2005 = 07/16/2005</li> <li>• dd/mm/yyyy example: 16 July 2005 = 16/07/2005</li> </ul> |
| Time Format       | Displays the time in one of the following styles: <ul style="list-style-type: none"> <li>• 12-hour example: 5:00 PM</li> <li>• 24-hour example: 17:00</li> </ul>                                              |



# 4

## Spot Vital Signs LXi Technical Overview

---

### System Description

The Spot Vital Signs LXi automatically measures systolic and diastolic pressure (excluding neonates), pulse rate, temperature (oral, adult axillary, pediatric axillary, rectal, and ear), and pulse oximetry (SpO<sub>2</sub>) as well as calculates Mean Arterial Pressure (MAP). Furthermore, Spot Vital Signs LXi allows the entry of height, weight, respiration rate, and pain level. Spot Vital Signs LXi also calculates Body Mass Index (BMI) following height and weight entry.

The device is intended to be used by clinicians and medically qualified personnel. It is available for sale only upon the order of a physician or licensed health care provider.

Refer to the Spot Vital Signs LXi Directions for Use manual for complete information.

### Battery System

Spot Vital Signs LXi uses a fuse and a shunt diode on the main board as well as a fuse in the battery assembly for protection against reverse battery connections. The fuse is located on the main board in the negative battery line. In the event of a reverse battery connection, a blown fuse will require replacement. The negative line location is more likely to short the battery (+) wire to ground than short the battery (-) wire to the VBUS node.

### Battery Charger

Provides two levels of power to fast-charge or to float the battery at fixed voltages using software control. The fast-charge level rapidly brings the battery up to near-full capacity (7 Volts) using software to not overcharge the battery. The float level, set per the manufacturer's recommendation, slowly tops-off the battery and can stay on indefinitely.

### "CPU I/O" Power Supply

Provides power to the I/O pins of the DragonBall MPU and the attached peripheral IC's. The nominal voltage is set to a value of 3.15 Volts so that both the 3.0 Volt MPU and 3.3 Volts peripherals can operate within specifications. Spot Vital Signs LXi uses buck switching topology for efficiency and an on-off power latch directly controls the power supply.

## Standby Mode

The Standby Mode conserves battery power. The device goes into Standby Mode if it is not used for two minutes. Press any button to bring the Spot Vital Signs LXi out of Standby Mode.

## CPU Core Power Supply

Powers the MPU core only and operates at the recommended Dragonball voltage range of 1.8 +/- 0.1 Volts; 2.0 Volts absolute max.

The +3.15 Volt supply powers the +1.8 Volt supply for efficiency and, more importantly, for power sequencing. During power-on the +1.8 Volt supply cannot exceed the +3.15 Volt supply. The R207/C182 further delays the +1.8 Volt supply. During power-off, D106 prevents the +1.8 Volt supply from exceeding the +3.15 Volt supply by more than a partial diode drop.

## +5V Supply

Provides supply to miscellaneous loads that require more than 3.15 Volts. The Power Latch controls the supply and the average load current is small, so efficiency is not important.

## Clock/Calendar Power

The Clock/Calendar U19 receives power from two sources. During power on the clock runs from the +5 Volt supply and during power down the clock receives power directly from the battery potential using a voltage divider. Capacitor C131 provides power to the U19 during the battery exchanges. Zener D102 prevents overvoltage to U19 in the case of operation directly from the charger.

## Mod F NIBP Power

This module is powered directly from VBUS through an EMI filter. Software via BP-PWR-ON-BAR controls the power. Level shifter Q115 / Q104 assures that the module remains off when the MPU is unpowered.

## SpO<sub>2</sub> and Thermometer Power

A common structure supplies SpO<sub>2</sub> and thermometry. A single push-pull driver delivers square wave power to a pair of isolation transformers, one for SpO<sub>2</sub> and the other for thermometry. THERM-PWR-BAR activates the driver under software control. Level shifter Q118, Q119 assures that the module remains off when the MPU is unpowered.

In the push-pull (forward) configuration, VBUS is multiplied by the transformer turns ratio and supplied to the input of the rectifier diodes. This is filtered and passed through low-noise 5-Volt linear regulators to supply power to the parameter modules.

## LCD Power

The LCD data interface connects directly to the +3.15 Volt MPU I/O power supply. This powers the driver ICs inside the LCD pane and assures the data interface to the LCD is always powered when the unit is powered.

The contrast (bias) supply to the LCD panel and the LCD CCFL backlight both derive power indirectly from a pre-regulator. A software control can shut down the pre-regulator. For efficiency, it is implemented as a buck switcher operating directly from the battery. The primary purpose of the pre-regulator is to protect the backlight inverter / ballast from battery voltages that it cannot tolerate and to control current out of the unregulated ballast. The voltage is set to the low side of 5 Volts to reduce battery consumption, while still driving the CCFL with adequate voltage.

## External Load (RS-232) Power

There are two power outlet means on the Spot Vital Signs LXi – RS-232 port power and the Braun PRO 4000 charger.

Two isolated power supplies provide power to each of the external RS-232 ports. This power is available to external devices on pin 9 of each serial connector whenever the software enables the corresponding RS-232 port.

## SureTemp Plus Thermometer Interface

The SureTemp Plus OEM module is located in an optional pod and is connected to the main PCB via a flex cable. The Thermometer Probe connection is located on a separate PC board. It is connected to the SureTemp Plus OEM module by a ribbon cable. The electrical interface to the thermometer pod is a 6-pin ribbon connector, Molex 52271-0690 or equivalent.

## Braun ThermoScan PRO 4000 Thermometer Interface

The interface to the Braun ThermoScan PRO 4000 thermometer is located in the optional pod and connects to the main PCB via a flex cable. The communications link is a single duplex current loop that also charges the Braun battery pack.

The presence or absence of current passing through the Braun ThermoScan PRO 4000 charge contacts controls the interface charge LED, located within the pod, which mimics the Braun ThermoScan Pro 4000 Base Station (type 6021).

## SpO<sub>2</sub> Description

A shielded flat-flex cable connects the user interface sensor connector J8 to the SpO<sub>2</sub> electronics.

The Masimo SpO<sub>2</sub> module mates to the Masimo Adaptor Board that, in turn, mates to the Main board. This provides power and communications to the module from the Main board.


The Nellcor SpO<sub>2</sub> module mounts directly to the Main board. This provides power and communications to the module, and brings the sensor signals to the main board where EMI and ESD can be controlled. The sensor signals pass through the main board to the mating connector for the shielded flat-flex cable.

# Radio Description

For the wireless radio description, refer to the manufacturer’s Directions for Use.

The radio interfaces to the Spot Vital signs LXi through an RS-232 female connector with a 9-pin male/female null modem connection.

**Figure 8. 9-pin Null Modem Pinout**



MA 0051 (DB9M/F) Adapter Pin out:

| 09 Male | 09 Female |
|---------|-----------|
| 2       | 3         |
| 3       | 2         |
| 4       | 1 + 6     |
| 5       | 5         |
| 1 + 6   | 4         |
| 7       | 8         |
| 8       | 7         |
| 9       | 9         |
| G       | G         |

# Bar Code Scanner Description

The Spot Vital Signs LXi has a 2D scanner (IMAGETEAM™ 4600). This scanner is programmed specifically for Welch Allyn. Scanners not purchased through Welch Allyn will not work with the Spot Vital Signs LXi.

## Interface

The bar code scanner interfaces to Spot Vital Signs LXi through an isolated RS-232 DB-9 female connector. Use Welch Allyn part number 4500-910 for the barcode option.



# Printer Description

The SPOT Vital Signs LXi interfaces with printer kit 4500-500.

**Table 10. LED Pattern Table**

| Pattern                                                         | Battery               | Buffer Mode                 |
|-----------------------------------------------------------------|-----------------------|-----------------------------|
| Constant Green                                                  | Running               | Normal                      |
| Fast Flashing Green                                             | Running               | Spool *                     |
| Long Green - Short Red                                          | Battery Low           | Normal                      |
| Short Green - Short Red                                         | Battery Low           | Spool *                     |
| Short Green - Long Orange                                       | Trickle Charging      | Normal                      |
| Long Orange Flash                                               | Trickle Charging      | Spool *                     |
| Fast Flashing Orange                                            | Fast charging         | Spool (Printing prohibited) |
| Fast Flashing Red                                               | Error Condition       | Spool (Printing prohibited) |
| No Light                                                        | Flat or in sleep mode |                             |
| * Spool may have resulted from Paper Out or Head Up conditions) |                       |                             |

No hardware power switch is fitted, as power control is either automatic or by command from the host computer. To save power, the printer enters sleep mode after a period of inactivity, factory-set to 30 seconds.

The printer external power supply is used to charge the printer batteries.

# Interconnect Diagram

(For drawing details see “Interconnect Diagram 1” on page 147 and see “Interconnect Diagram 2” on page 148.

| Item   | Description                                   | Item  | Description                                 |
|--------|-----------------------------------------------|-------|---------------------------------------------|
| A1     | PCA, Main, Spot Vital Signs LXi               | W1    | Assy, Power Harness, Spot Vital Signs LXi   |
| A1-F1  | Fuse, FlastBlo, 5A                            | W1-A1 | PCA, Power Interface, Spot Vital Signs LXi  |
| A2     | PCA, Serial Interface, Spot Vital Signs LXi   | W2    | Flex, Main, Spot Vital Signs LXi            |
| A3     | Display, LCD, 320 x 240, Mono                 | W3    | Cbl Flx, 1/2mm, 15 POS                      |
| A4     | Inverter, LCD                                 | W4    | Assy, Ballast Harness, Spot Vital Signs LXi |
| A5     | Assy, Battery, Spot Vital Signs LXi           | W5    | Assy, LCD Gnd Harness                       |
| A5-F1  | Wire Protection Fuse                          | W6    | Assy, Cable, Flat Flex, SpO <sub>2</sub>    |
| A6     | Assy, Fan Unit Cooling                        |       | Wire Clip, Side Entry,.375 Wide             |
| A7     | PCA, Masimo Adapter, Spot Vital Signs LXi     |       | Foam Strip, 2 Sided Adhesive                |
| A8     | Masimo MS-11 SpO <sub>2</sub>                 | W7    | Cbl, Flx, 1mm, 6POS                         |
|        | Shield, Masimo EMI                            | W8    | Cbl, Flx, 1mm, 8POS                         |
|        | Tape, Copper, Masimo Shield                   | W9    | Cbl, Flx, Braun, Spot Vital Signs LXi       |
| A9     | PCA, SureTemp Plus, OEM w/o EDAC              |       |                                             |
| A10    | PCBA, BP OEM, (MOD F)                         |       |                                             |
| A11    | Assy, Pump, and Tubing, Spot Vital Signs LXi  |       |                                             |
| A11-F1 | Filter, Air, 1/8 Barbs                        |       |                                             |
| A12    | Assy, Transducer Tubing, Spot Vital Signs LXi |       |                                             |
| A13    | PCA, EDAC, Spot Vital Signs LXi               |       |                                             |
| A14    | Nellcor MP506 SpO <sub>2</sub> PCB            |       |                                             |
| A15    | PCA, Braun Interface, Spot Vital Signs LXi    |       |                                             |

## Wire Color

|     |        |
|-----|--------|
| BLK | Black  |
| BRN | Brown  |
| RED | Red    |
| ORN | Orange |
| GRN | Green  |
| BLU | Blue   |
| PNK | Pink   |

# 5

## Maintenance and Service

---

### Cleaning



**WARNING** Before performing any maintenance or service to the Spot Vital Signs LXi, disconnect the AC power line from the electrical outlet.

### Spot Vital Signs LXi



**Caution** Prevent water or other fluids from entering any connectors. If the connectors get wet, dry them with warm air. Check all measurement functions.

**Caution** Do not sterilize or autoclave the Spot Vital Signs LXi device.

As necessary, clean the device with a cloth slightly dampened with either 70% isopropyl alcohol, 10% chlorine bleach solution, or mild detergent in water, or use PDI Sani-System Cloths. Never immerse Spot Vital Signs LXi in any type of fluid.

### Blood Pressure Cuff

Before washing the Reusable Two-Piece Cuff, remove the two connectors and bladder, close off tubes with plugs (available as accessory 5082-163), and close the hook and loop fasteners. After washing, allow the cuff to air dry. Re-assemble the tube fittings.

Never press with an iron.

### Blood Pressure Hose and Cable

Wipe the pressure hose with a damp cloth moistened in a mild detergent solution. Do not immerse hose.

## SureTemp Plus Thermometer

### Temperature Probe



**Caution** DO NOT immerse or soak the probe in any type of fluid.

**Caution** DO NOT use steam, heat, or gas sterilization on the probe.

**Caution** DO NOT autoclave the probe.

Press down on the connector tab and slide the connector out of the port to remove the temperature probe.

Regularly wipe the probe with a cloth dampened with warm water and a mild detergent solution, a 70% isopropyl alcohol solution, or a 10% chlorine bleach solution.

### Removable Probe Well



**Caution** DO NOT use hard or sharp objects to clean the probe well. This could damage the probe well and cause the device to not function properly.

**Caution** DO NOT use steam, heat, or gas sterilization on the probe well.

**Caution** DO NOT autoclave the probe well.

1. Remove the temperature probe from Spot Vital Signs LXi (see “Temperature Probe” on page 44).
2. Grasp the well under the probe opening and pull up gently to remove it from the device.
3. Swab the inner and outer surface of the probe well with a cloth dampened with a mild detergent solution, 70% isopropyl alcohol, or 10% chlorine bleach solution. Immerse the probe well in mild detergent solution as necessary for cleaning.
4. Dry all surfaces thoroughly before re-assembling the device (see “SureTemp Plus” on page 26) for reassembly instructions.

## Braun ThermoScan PRO 4000 Thermometer

Use a soft cloth slightly moistened with alcohol to clean the thermometer display and exterior. Do not use abrasive cleaners.

Damage to the probe tip or the presence of dirt or cerumen on the probe tip can affect the accuracy of your temperature measurement. To clean the tip, gently wipe it with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. Allow to dry at least five minutes before taking a temperature.

Every month, clean the Braun ThermoScan PRO 4000 charging contacts within the Spot Vital Signs LXi housing with a swab slightly dampened with alcohol.

## SpO<sub>2</sub> Sensors



**WARNING** Do not immerse the sensor in water, solvents, or cleaning solutions (the sensors and connections are not waterproof). Do not use irradiation, steam, or ethylene oxide for sterilization.

Clean the reusable SpO<sub>2</sub> sensor with a 70% isopropyl alcohol solution. Do not immerse the sensor.

## Radio

Refer to Manufacturer's Directions for Use.

## Printer

The printer label kit comes with a two-step head cleaner. Follow the directions provided with the cleaner.

# Battery Replacement

## Spot Vital Signs LXi




**Caution** Only use the Welch Allyn 4500-84 lead acid battery. Using an incorrect battery will cause damage to the Spot Vital Signs LXi and void the warranty.

**Caution** Do not break the shrink-wrap around the battery.

If necessary, replace the internal battery after heavy use or the battery no longer charges.

1. Power off the Spot Vital Signs LXi and disconnect the AC power transformer cord.
2. Remove the four screws holding the battery door using a phillips-head screwdriver. Remove the battery door to expose the battery.
3. Tip the Spot Vital Signs LXi backward and slide the battery out. Disconnect the one-way connector and then attach a new battery to the connector as shown. The one-way connector ends only connect one way. Do this as quickly as possible to prevent loss of clock time.

**Figure 9. Battery Removal and One-way Connector**



4. Slide the battery into the compartment as far as it will go. Push the connector down into the case next to the battery.
5. Replace the battery door and tighten each of the four screws.
6. Connect the AC power transformer to the Spot Vital Signs LXi and allow the new battery to charge for approximately 9 hours. The rechargeable batteries in the Braun ThermoScan PRO 4000 thermometer require an additional 1 hour to charge. You can use the Spot Vital Signs LXi during this charging period via the AC power cord.

If Spot Vital Signs LXi displays the E38 error code after power up, set the date (see "Date/Time Menu Options" on page 33).



The battery is a non-spillable lead-acid battery and must be recycled or disposed of properly according to national or local regulations.

## Braun ThermoScan PRO 4000



**Caution** Do not use alkaline batteries in the Braun ThermoScan PRO 4000 thermometer.

Welch Allyn supplies a rechargeable battery pack with the Braun ThermoScan PRO 4000 thermometer.

1. Open the battery compartment.
2. Remove the battery pack and replace with a new battery pack, verify the poles are in the right direction.
3. Slide the battery door back in until it snaps into place.

If the battery is completely discharged, the LEDs will not illuminate. Allow the unit to charge at least 15 minutes before proceeding.



The battery is a rechargeable battery and must be recycled or disposed of properly according to national or local regulations.

## Printer



**Caution** Do not use alkaline batteries in the printer.

1. Open the battery compartment.
2. Remove the battery pack and replace with a new battery pack, verify the poles are in the right direction.
3. Slide the battery door back in until it snaps into place.

Allow the printer to charge at least 9 hours before proceeding.



The battery is a rechargeable battery and must be recycled or disposed of properly according to national or local regulations.

## Spot Vital Signs LXi Firmware Upgrade CD

This CD upgrades the Spot Vital Signs LXi unit software version. A detailed procedure is available on the upgrade CD.

## Spot Vital Signs LXi Radio Configuration CD

The software on this CD configures the Spot Vital Signs LXi Radio for use on a local wireless network.

## Product Disposal

Do not dispose of this product as unsorted municipal waste. Prepare this product for reuse or separate collection as specified by Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste Electronic and Electrical Equipment (WEEE). If this product is contaminated, this directive does not apply.

For more specific disposal information, see [www.welchallyn.com/weee](http://www.welchallyn.com/weee), or contact Welch Allyn Customer Service at +44 207 365 6780.



# 6

## Calibration and Verification

---

### Terms

**BP test volume repair fixture.** Equipment with four separate 100cc, 250cc, and two 500cc linearity volume cylinders. Referred to in this procedure as “Test Volume”.

**Calibration.** Process to change the outputs of the UUT that are out of range. After saving, the changes become permanent until the next calibration.

**Pneumatic test assembly.** Equipment that connects the UUT, the pressure meter, the squeeze bulb, and the BP test volume repair fixture.

**Spot LXi Service Tool.** Software package to verify and calibrate the Spot LXi device.

**Service computer.** Computer with the repair software loaded and ready to operate.

**Functional Tests.** Process to check the outputs of the UUT.



**UUT.** Unit under test.

## Required Test Equipment

| Description                                                                                                                                | Part Number | Qty |
|--------------------------------------------------------------------------------------------------------------------------------------------|-------------|-----|
| Power supply, 0-10 Volts, 2 amps capacity or better, fine voltage adjustment required                                                      |             |     |
| Digital multi meters (DMM)                                                                                                                 |             | 2   |
| Battery simulator connector, 24 inch or shorter                                                                                            |             |     |
| Patch cords- banana plug to banana plug, Red, 18 AWG, 24 inch or shorter                                                                   |             | 3   |
| Patch cords- banana plug to banana plug, Black, 18 AWG, 24 inch or shorter                                                                 |             | 2   |
| USB Cable, 2.0/5-Pin Mini-B Cable                                                                                                          |             |     |
| Calibrated pressure meter                                                                                                                  |             |     |
| PC running Spot Lxi Service Tool (Ver 3.0 or greater)                                                                                      | 4500-905    |     |
| 9600 Plus Temperature Calibration Tester                                                                                                   | 01802-110   |     |
| Nellcor-approved SpO2 simulator, <i>for devices with Nellcor option</i>                                                                    | SRC-MAX     |     |
| Masimo SpO2 Tester, <i>for devices with Masimo option</i>                                                                                  | 1795        |     |
| Thermometer Calibration Key, Optional                                                                                                      | 06138-000   |     |
| Blood pressure two-tube hose with Y fitting                                                                                                | 5082-183    |     |
| Blood pressure two-tube hose (standard cuff)                                                                                               | 4500-30     |     |
| BP test volume repair fixture (includes items below)<br><b>Note</b> This is the new Test Volume needed for passing the NIBP Linearity Test | 407672      |     |
| BP test volume repair fixture. Must have at least the following volume cylinders: 100cc, 250cc, and 500cc                                  |             |     |
| Squeeze bulb with bleeding screw                                                                                                           |             |     |
| Square (4-way) tube fitting connector                                                                                                      |             |     |
| 12-inch (30.5 cm) tubing with .25-inch (6.35 mm) ID                                                                                        |             | 4   |

## Pneumatic test assembly setup


1. Connect the squeeze bulb with bleed screw to a square (4-way) tube fitting connector, a calibrated pressure meter, a test volume, and to a blood pressure two-tube, Y-fitting hose as shown in [Figure 1](#) on page 51.

**Figure 1. Pneumatic test assembly setup**

## Power Supply Test Setup


1. Connect the power supply to the Battery Simulator connector, banana plug patch cords, digital multi meter (DMM) for voltage monitoring, DMM for current monitoring.

**Figure 2. Power Supply Test Setup**



Test the power supply output as shown in [Figure 3](#). Observe the correct polarity.

**Figure 3. Power Supply**



2. Turn on the power supply and adjust to 6.5 +/- 0.1V DC before connecting the UUT. Turn off power supply.
3. Disconnect the UUT battery and connect the external power supply to the UUT, using the Battery Simulator. Observe the correct polarity.



**Caution** DO NOT connect the Spot LXi AC power adapter to the UUT power connection port. This will result in damage to the UUT, AC power adapter, or the Battery Simulator power supply.

4. Turn on the Power Supply.

**Note** All tests and verification in this procedure require battery voltage input adjusted to 6.5V DC; Except as called out in the Voltage Calibration section.

## Battery Removal / Power Shutdown

Removal of the UUT battery, or power shutdown of the Battery Simulator external power supply, will require reset of the date and time prior to any further testing, calibration, or verification activities.

### Time & Date Reset:

1. Power on the UUT in configuration mode by pressing the **Check** and **Power On** buttons on the UUT. Hold the **Check** button until the UUT beeps twice. The UUT screen displays any POST errors detected.

**Figure 4. POST Errors Detected**



2. Click the **Check** button on the UUT. The Configuration Menu opens.
3. Scroll to **Date/Time** on the UUT Configuration Menu and select.

**Figure 5. Date/Time Menu**

4. Scroll to **Date** on the UUT and select. Using the navigation keys, adjust the correct date and select. Scroll to **Time** and select. Using navigation keys, adjust the correct time and select.

**Figure 6. Set Date, Set Time**

5. Scroll to **Exit** and select.

**Figure 7. Exit Menu**


6. Reset the UUT using the power button. Do NOT turn off the power supply. The next power up of the UUT screen displays the correct date and time without errors.

## Voltage Calibration

1. Remove the battery from the UUT. Connect the UUT to a Battery Simulator with the external power supply set to 6.5 +/- 0.1V DC. (Follow the steps in "[Power Supply Test Setup](#)" on page 52.)
2. Connect the service PC running the Spot Lxi Service Tool to the UUT with the USB Cable. Open the Spot LXi Service tool from the PC.
3. Power on the UUT in configuration mode by pressing the **Check** and **Power On** buttons on the UUT. Hold the **Check** button until the UUT beeps twice.

**Figure 8. Configuration Menu**

4. Verify the power supply is at  $6.5 \pm 0.1$  V DC at the battery connector simulator leads.
5. Verify that the UUT is connected and communicating with Spot Lxi Service Tool.  
**Select Unit > Information.** The Unit information window opens. Click **OK** to close.
6. Select **Calibrate > Voltage** on the service computer. The UUT screen goes blank.

**Figure 9. Calibrate Voltage**


7. Carefully adjust the power supply to  $5.5 \pm 0.1$  V DC at the battery connector simulator leads. A fine voltage adjustment will be needed to make this adjustment.



**Note** Do not adjust power supply below 5.5V DC or the UUT will go into shutdown mode.

8. Read the voltage on the DMM that is connected to the battery simulator leads.
9. Enter this voltage reading into the **Actual Voltage** field and click the **Calibration** button.
10. Click the **OK** button to accept the inputs. Click **Yes** to make the changes permanent.
11. Enter your initials in the Calibration Signature field and click **OK**. The UUT resets and powers up in normal mode.

**Figure 10. Calibration Signature Field**



**Note** After "Reset" and power up, the UUT detects the low 5.5V DC supply setting and goes into auto shutdown. Adjust power back to 6.5V DC.


## Pressure verification test

Spot LXi devices require annual calibration.

**Note** If this is a new NIBP board, go to "[NIBP board Initialization](#)" on page 62 and carefully follow the instructions to initialize this board.

1. Connect the UUT to the pneumatic test assembly using the two-tube blood pressure hose with the Y-fitting.

**Figure 11. Blood Pressure Hose with Y-fitting**



2. Connect the test volume hose to the 500cc cylinder.

3. Power on the UUT in configuration mode by pressing the **Check** and **Power On** buttons on the UUT. Hold the **Check** button until the UUT beeps twice. The configuration screen displays.

**Figure 12. Configuration Menu**



Scroll to Blood Pressure and select it


**Figure 13. Blood Pressure Menu**




Select **BP Calibration Check**.

**Figure 14. BP Calibration Check**

4. Open the Spot LXi Service Tool on the service computer and select **Test > Calibration**. A Test Calibration box displays on the screen.

**Figure 15. Test Calibration with Open Valve**

5. Verify the pressure is at or near zero and the voltage is within 0.2V DC of the power supply (6.5V DC). Verify the valve status is open and the pump status is off.
6. Select **Close Valve** on the UUT. The valve status changes to closed.

**Figure 16. BP Calibration Check Close Valve****Figure 17. Test Calibration with Closed Valve**

7. Select **Start Cuff Inflation**. The pump starts inflating the cuff. Verify the pump status is on and the pressure is rising.

**Figure 18. Start Cuff Inflation**

8. Inflate the pressure up to 250 mmHg. Select the **Stop Pump** on the UUT to stop the inflation. Fine tune the pressure reading to 250 mmHg with the squeeze bulb.

**Figure 19. Fine Tune the Pressure Reading to 250 mmHg**

9. Check the reading on the pressure meter against reading on the UUT.
  - a. Verify the pressure on the UUT is within the pressure tolerance in the table. Use the squeeze bulb to fine tune the pressure at each inflation pressure limit. If the pressure is within tolerance, go to the next inflation pressure.

| Inflation pressure | Pressure tolerance | Inflation pressure | Pressure tolerance |
|--------------------|--------------------|--------------------|--------------------|
| 250 mmHg           | +/- 2.0 mmHg       | 50 mmHg            | +/- 1.0 mmHg       |
| 150 mmHg           | +/- 1.5 mmHg       | 0 mmHg             | +/- 1.0 mmHg       |

- b. If all readings are within specification, go to [“Functional tests” on page 67](#).
- c. If the UUT is not within the pressure tolerance, open the valve to bleed pressure to zero. Select **Calibrate > Pressure** on the repair software. Proceed to [“Blood pressure calibration” on page 63](#). UUT needs calibration.
- d. If replacement of the NIBP board is required, first proceed to [“NIBP board Initialization” on page 62](#) and then proceed to [“Blood pressure calibration” on page 63](#).

## Blood pressure calibration

**Note** Complete the voltage calibration before beginning a blood pressure calibration.

If all readings were within specifications during [“Pressure verification test” on page 57](#), skip this section and go to [“Functional tests” on page 67](#).

## NIBP board Initialization

**Note** All New NIBP Boards must be initialized before calibration.

1. Power on the UUT in configuration mode by pressing the **Check** and **Power On** buttons on the UUT. Hold the **Check** button until the UUT beeps twice. The configuration screen displays.


**Figure 20. Configuration Menu**





2. Open the Spot LXi Service Tool on the service computer and select **Calibrate > Reset Defaults**. The Reset Defaults window displays.

**Figure 21. Reset Defaults Window**



3. Click **Reset > Initialize > OK**. The NIBP board is initialized.
4. Perform the "Voltage Calibration" on page 55 before proceeding to "Blood pressure calibration".

### Blood pressure calibration

1. Power on the UUT in configuration mode by pressing the **Check** and **Power On** buttons on the UUT. Hold the **Check** button until the UUT beeps twice.

**Figure 22. Blood Pressure Menu**




2. The configuration screen displays, scroll to **Blood Pressure** and select it. Select **BP Calibration Check**.

**Figure 23. BP Calibration Check**




3. Open the Spot LXi Service Tool on the service computer and select **Calibrate > Pressure**. The Calibrate Pressure window opens.


**Figure 24. Calibrate Pressure Window**






4. Connect the UUT to the pneumatic test assembly using the two-tube blood pressure hose with the Y-fitting. (See “[Blood Pressure Hose with Y-fitting](#)” on page 57.)
5. Connect the pneumatic test assembly hose to the 500cc cylinder of the test volume and verify the following:
  - the pressure is at zero mmHg,
  - the bleed screw on squeeze bulb is closed,
  - valve on UUT is in the open state.
6. On the Spot LXi Service Tool select **Calibrate Zero button**. 
7. Close Valve on the UUT and use the squeeze bulb to manually inflate the UUT to 200 mmHg. Let stabilize.

**Figure 25. Fine Tune the Pressure Reading to 200 mmHg**




8. Enter the pressure meter reading in the **Calibration Gain** field on the service computer. Be as precise as possible. If pressure meter has resolution of 1.00 mmHg make sure it reads 200.00 without switching to 199.00 or 201.00. Enter 200.00 in **Calibration Gain** field.

**Figure 26. Calibration Gain Field**

The 'Calibrate Pressure' dialog box has a blue title bar with a close button. It contains a 'Status:' section with three read-only fields: 'Internal Temperature: 80.2 °F', 'Pressure: 0.00 mmHg', and 'Valve Status: Open'. Below this is a text instruction: 'Ensure that valve is open, and that no residual pressure is applied to unit, before pressing "Calibrate Zero" button.' followed by a 'Calibrate Zero' button. Another instruction follows: 'Enter calibrated pressure reading (from external meter) in space below, then press "Calibrate Gain" button.' followed by an empty text input field and a 'Calibrate Gain' button. At the bottom are three buttons: 'OK' with a green checkmark, 'Cancel' with a red X, and 'Help' with a question mark.

9. Click **Calibrate Gain** on the service computer.
10. Click **OK** to accept the inputs. Click **Yes** to make the changes permanent.
11. Enter your initials in the **Calibration Signature** field and click **OK**. The UUT resets and powers up in normal mode.

**Figure 27. Calibration Signature Field**

The 'Calibration Signature' dialog box has a blue title bar with standard window controls. It contains the instruction: 'Please enter a three character calibration signature, eg. Your initials DED'. Below this is a text label 'Calibration Signature:' followed by an empty text input field. At the bottom are two buttons: 'OK' with a green checkmark and 'Cancel' with a red X.

12. Blood Pressure calibration is complete.

# Functional tests

**Note** Perform the Functional Tests sequentially. Fill out the “[Service Work Checklist](#)” on page 111 with the Calibration and Functional Test data obtained as you perform these steps.

Each test is software driven. Selecting the Test button begins the test.

Connect the correct Test Volume cylinder called out in the procedure (in parenthesis below) or the test will fail.

## Test specification details:

|                            |                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Leak test</b>           | 6mm maximum in 15 seconds with a volume of <b>(100cc)</b> at 250 mmHg                                                                                                                                                                                                                                                                                        |
| <b>Dump test</b>           | From 260 to less than 15 mmHg in 10 seconds or less with a volume of <b>(500 cc)</b>                                                                                                                                                                                                                                                                         |
| <b>Inflation test</b>      | From 5 to 210 mmHg in 7 seconds or less with a volume of <b>(250 cc)</b>                                                                                                                                                                                                                                                                                     |
| <b>Valve control test</b>  | From 160 mmHg with a volume of <b>(100 cc)</b><br>Open the valve for the following time intervals and pressure drops: <ul style="list-style-type: none"> <li>• 10 msec = 4 to 12 mmHg</li> <li>• 15 msec = 4 to 15 mmHg</li> <li>• 25 msec = 4 to 25 mmHg</li> </ul>                                                                                         |
| <b>Inflation linearity</b> | Connect a standard cuff hose between the UUT and the 500cc Linearity input on the Test Volume and inflate. <ul style="list-style-type: none"> <li>• Measure the pressure at 6, 9, and 12 seconds and compare to the expected readings at these intervals.</li> <li>• Also check the noise level at 7, 10, and 13 seconds into the inflation step.</li> </ul> |


## Functional test steps:

1. If the UUT is connected to the Battery Simulator, set the external power supply to 6.5 +/- 0.1V DC.

**Note** The Spot LXi AC power adapter should NOT be connected to the UUT.

If the UUT is NOT connected to the battery simulator, connect the Spot LXi AC power adapter to the UUT.



2. Power on UUT in normal mode.
3. Connect the UUT to the pneumatic test assembly using the two-tube blood pressure hose with the Y-fitting. (See “[Blood Pressure Hose with Y-fitting](#)” on page 57.) Verify that the bleed valve on the squeeze ball is closed.
4. Open the Spot LXi Service Tool on the service computer and select **Test > Pneumatics**. The Test Pneumatics window opens.

**Figure 28. Test Pneumatics Window**


- a. Connect the test specific cylinder in parenthesis; (100cc), (500cc), (250cc)
- b. Sequentially select each pneumatic test, starting with "Leak (100cc)". Click **Test**. A successful test will display "Passed".

**Figure 29. Passed Leak Test**

- c. Perform each of the tests listed after connecting the test specific cylinder. The "Inflation Linearity (Cuff)" test may use the 500cc cylinder in place of the adult cuff.

**Figure 30. Inflation Linearity Cuff Test****Figure 31. Inflation Linearity 500cc Cylinder**


- d. Click **OK** after completing the last test to close the **Pneumatics** test window.
5. On the Spot LXi Service Tool select **Test > Current Levels**. Perform the tests listed. The Test Current Levels window opens.

**Figure 32. Test Current Levels Window**


- a. **Sleep State > Test**, screen is off and the UUT is idle. Verify that the current is less than 500 mA. Click **OK** after completing the measurement.

**Figure 33. Sleep State**

- b. **Valve/pump Mode > Test**, the UUT is fully active with the valve closed and the pump running. Verify that the current is less than 2.0 amperes in this state. Click **OK** after completing the measurement.
- c. Normal Mode; screen is lit, but the UUT is inactive. The UUT current draw is less than 1.25 amperes in this state.
- d. Click **OK** after completing the last test to close the **Test Current Levels** window.
6. On the Spot LXi Service Tool select **Test > Noise Levels > Test**. The **Test Noise Levels** window opens.

**Figure 34. Test Noise Levels Window**

- a. The Test Noise Levels test measures the amount of noise on the UUT A-D pressure channel. The sample is 1 second in length with 0 mmHg applied. The maximum allowed noise is .050 mmHg.
  - b. Click **OK** after completing the test to close the **Test Noise Levels** window.
7. On the Spot LXi Service Tool select **Test > Buttons**. The **Test Buttons** window opens.

**Figure 35. Test Buttons**


- a. Press each button on the UUT and verify the appropriate button lights on the Spot LXi Service Tool window.

Figure 36. UUT Buttons



- b. Click **OK** after completing the test to close the **Buttons Test** window.
- 8. On the Spot LXi Service Tool select **Test > Interface**. The **Test Interface** window opens.

Figure 37. Test Interface Window



- a. Perform all commands in this test, verifying the screen on the UUT responds correctly with each command:

|                                                   |  |
|---------------------------------------------------|--|
| Normal                                            |  |
| All On                                            |  |
| All Off                                           |  |
| Test Pattern                                      |  |
| Backlight On/Off, = blank dark UUT display screen |  |



|                                                        |  |
|--------------------------------------------------------|--|
| Power Indicator On/Off                                 |  |
| Charge Indicator On/Off                                |  |
| Enunciator On/Off, Toggles enunciator tone on and off. |  |

Figure 38. Normal



Figure 39. All On



**Figure 40. All Off**



**Figure 41. Test Pattern**



**Figure 42. Backlight On/Off****Figure 43. Power Indicator On/Off****Figure 44. Charge Indicator On/Off**

- b. Click **OK** after completing the tests to close the **Test Interface** window.

**To perform the over pressure test:**

1. Verify that the UUT is in normal mode.
2. Re-install the pneumatic test assembly with the 250cc cylinder connected.
3. On the Spot LXi Service Tool select **Unit >Status**. A **Unit Status** window displays.

**Figure 45. Unit Status Window**

The screenshot shows the 'Unit Status' window with two tabs: 'Status Information' and 'POST Errors'. The 'Status Information' tab is active. It contains the following sections:

- Unit Status:**
  - Powerup Mode: **Configuration**
  - Charge State: **Battery**
  - Battery Level: **Above 6.15 Volts**
- ModF Status:**
  - Module State: **Dumping**
- Hardware/Patient States:**
  - Valve Status: **Open**
  - Pump Status: **Off**
  - Malfunction (E11): **No**
- Safety States:**
  - Over Pressure: **No**
  - Autozero Failure: **Yes**
  - Battery Failure: **No**
  - Ambient Out Of Range: **No**
  - Unrecoverable POST: **No**
  - Residual Pressure Failure: **No**
- Cuff Pressure:** **0.00 mmHg**
- ☐ **Log Status Information** (This checkbox is circled in blue in the original image)

At the bottom of the window are two buttons: **OK** (with a green checkmark icon) and **Help** (with a question mark icon).

4. Check the **Log Status Information** box in the lower right hand corner of the Spot LXi Service Tool. The Event logs display in rapid succession.
5. Immediately initiate a blood pressure reading. As the pressure increases, start applying additional pressure with the squeeze bulb.
6. Increase the pressure to the UUT until an E10 error displays.
7. Uncheck the **Log Status Information** box and click **OK** after completing the test to close the **Unit Status** box.

Scroll through the event log to locate the last 'no' answer, directly followed by the first 'yes' answer in the column labeled over pressure.

Figure 46. Event Log



| Battery Level | State     | Valve State | Pump State | E11 | Over Pressure | Auto Fail |
|---------------|-----------|-------------|------------|-----|---------------|-----------|
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | On         | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | Off        | No  | No            | No        |
| ATT_LEVEL_1   | Inflating | Closed      | Off        | No  | No            | No        |
| ATT_LEVEL_1   | Dumping   | Open        | Off        | No  | Yes           | No        |
| ATT_WARN      | Dumping   | Open        | Off        | No  | Yes           | No        |
| ATT_WARN      | Dumping   | Open        | Off        | No  | Yes           | No        |
| ATT_WARN      | Dumping   | Open        | Off        | No  | Yes           | No        |

8. Scroll to the right to see the mmHg pressure that these events occurred.
9. Select the higher pressure reading when the E10 error occurred. Verify the pressure is between 296 and 329 mmHg.

**Note** The following tests are for verification only. There is no updating of calibration parameters.

## Temperature Functional Check

### Temperature test setup


The 9600 Plus Calibration Tester takes approximately 20 minutes to heat to the lowest setting. When testing several thermometers at all three temperatures, it is recommended to test all probes at one Calibration Set Point Temperature before proceeding to the next Calibration Set Point Temperature.

To further expedite testing start at the lowest Calibration Set Point Temperature. The 9600 Plus Calibration Tester does not have an internal fan, this causes a longer cool down time than warm up time.

Refer to the 9600 Plus Calibration Tester Directions for Use manual for specific information regarding the LCD window or the control buttons.

1. Choose the proper mains plug insert and slide it over the two prongs in the power converter.

**Figure 47. Power Adapter and Mains Plug Inserts**



2. Plug the power adapter into the 9600 Plus Calibration Tester and the opposite end into a wall outlet.
3. Place the 9600 Plus Calibration Tester on a level surface away from sunlight, drafts, and other sources of heat or cold.
4. Observe the Set Point Mode in the upper left hand corner of the LCD display. If the unit displays a "D", it is in Default Mode and will heat to the lowest Set Point Temperature. If you do not want to conduct testing at this Set Point Temperature, press and hold the **Temperature Selection** button to select the desired setting. The temperature display will flash before staying on continuously to indicate the 9600 Plus Calibration Tester has stabilized and is ready for use.

## Changing the Calibration Set Point Temperature on the 9600 Plus Calibration Tester

To scroll from one set point to the next, press and hold the **Temperature Selection** button (up/down button) until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, starts to flash, and continues flashing until the cavity reaches the equilibrium at the new set point.

## Welch Allyn SureTemp Plus Thermometers




**Caution** Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow for thermal accommodation.

**To begin functional verification of the SureTemp Plus thermometer:**

1. Verify that the UUT is ON and is operating in normal mode.
2. Remove the thermometer probe from the probe well and clean it with either a 70% isopropyl alcohol solution, or a 10% chlorine bleach solution. Let the thermometer air dry. Do not apply a probe cover.
3. Verify that the blue (oral) probe well is in place.

Insert the thermometer into the 9600 Plus Calibration Tester device thermistor device port (small port).

**Figure 48. 9600 Plus Calibration Tester with the Welch Allyn SureTemp Plus Thermometer**



4. Verify that the thermometer enters Monitor Mode, (approximately one minute), the snail symbol and "Monitor" are displayed on the UUT display.

**Figure 49. Temperature Monitor Mode**



5. After the UUT enters "Monitor" mode wait for approximately one minute or until the temperature on the thermometer is stable for ten seconds.
  - a. Verify that the UUT thermometer's temperature reading is within the temperature tolerance in the table.

| Calibration Set Point Temperature | Temperature tolerance |
|-----------------------------------|-----------------------|
| 96.8°F (36.0°C)                   | +/- 0.2°F (+/- 0.1°C) |
| 105.8°F (41.0°C)                  | +/- 0.2°F (+/- 0.1°C) |



- b. When verifying several UUT thermometers at both temperatures, it is recommended to verify all UUT thermometers at one Calibration Set Point Temperature before proceeding to the next Calibration Set Point Temperature.

- c. Changing the Calibration Set Point Temperature (9600 Plus Calibration Tester):

To scroll from one set point to the next, press and hold the **Temperature Selection** button (up/down button) until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, starts to flash, and continues flashing until the cavity reaches the equilibrium at the new set point.

#### To begin functional verification of SureTemp Plus without thermometer (Optional):

If the UUT fails the temperature verification using the SureTemp Plus thermometer; then perform this separate functional test of the UUT temperature circuitry.

1. Remove probe and probe well completely.
2. Insert calibration key (part number 06138-000) and replace probe well.

**Figure 50. Calibration Key and Probe Well**



3. Insert and then remove probe from probe storage channel to reset thermometer.



**Figure 51. Reset Thermometer**

4. Wait for display test, and then observe display.
5. Display should read  $36.3^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$  or  $97.3^{\circ}\text{F} \pm 0.2^{\circ}\text{F}$ .

**Figure 52. Temperature Display in F°**

6. Remove calibration key and install probe well.
7. If the UUT temperature circuitry verified successfully and the temperature verification failed using thermometer probe; then replace SureTemp Plus thermometer probe.

## Braun ThermoScan Pro 4000



**Caution** Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow for thermal accommodation.

### To begin functional verification of the Braun ThermoScan PRO 4000 thermometer:

1. Clean the probe tip with a cotton swab slightly moistened with isopropyl alcohol, remove excess alcohol with a clean cotton swab, and let air dry for 5 minutes. Do not use any chemical other than alcohol to clean the probe window.
2. Place the Braun ThermoScan Pro 4000 in Calibration Mode using the following steps.
  - a. Push and release the **I/O mem** button to turn the product on. The display shows symbols and functions. The Braun ThermoScan Pro 4000 performs an automatic self check. After a few seconds you hear a beep and see two dashed lines, a sound symbol, and ° C or ° F on the display.

**Figure 53. Braun Self Check**



- b. Push the **I/O mem** button again and keep it pressed for the following steps:
  - After approx. 3 seconds the “OFF” symbol flashes on the display (keep pressing the button).
  - When you hear a beep, release the button immediately.
  - The Braun ThermoScan Pro 4000 is now in calibration check mode and the display is flashing and showing the “CAL” symbol.
  - Apply a new probe cover. Insert the Braun ThermoScan Pro 4000 into the 9600 Plus Calibration Tester ear device port (large port). Place the probe firmly into the ear device port.

**Figure 54. Braun Ear Device Port**



3. Wait approximately three seconds, press the thermometer **start** button, and watch for the “Exac Temp” light to flash.

**Figure 55. Press Start Button**



4. Leave the thermometer in the 9600 Plus Calibration Tester until you hear a beep. Remove the Pro 4000 Thermometer from the 9600 Plus Calibration Tester and read the temperature in the thermometer's display.
  - a. Verify that the UUT thermometer's temperature reading is within the temperature tolerance in the table.

| Calibration Set Point Temperature | Temperature tolerance |
|-----------------------------------|-----------------------|
| 96.8°F (36.0°C)                   | +/- 0.2°F (+/- 0.1°C) |
| 105.8°F (41.0°C)                  | +/- 0.2°F (+/- 0.1°C) |

- b. Wait one full minute before taking another reading with the same thermometer. Repeated measurements in short sequence may cause higher readings.
  - c. When verifying several UUT thermometers at both temperatures, it is recommended to verify all UUT thermometers at one Calibration Set Point Temperature before proceeding to the next Calibration Set Point Temperature.
  - d. Changing the Calibration Set Point Temperature (9600 Plus Calibration Tester):
 

To scroll from one set point to the next, press and hold the **Temperature Selection** button (up/down button) until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, starts to flash, and continues flashing until the cavity reaches the equilibrium at the new set point.
5. The Braun ThermoScan PRO 4000 thermometer will exit CAL mode after four minutes.

**To begin functional verification of the Braun ThermoScan PRO 4000 thermometer Down Load:**

1. Clean the probe tip with a cotton swab slightly moistened with isopropyl alcohol, remove excess alcohol with a clean cotton swab, and let air dry for 5 minutes. Do not use any chemical other than alcohol to clean the probe window.
2. Insert the Braun ThermoScan PRO 4000 thermometer into the UUT and wait for the "Exac Temp" light to go out.

**Figure 56. Braun Exac Temp Light**



3. Remove the Braun ThermoScan PRO 4000 thermometer from the UUT.
4. Apply a new probe cover. Insert the Braun ThermoScan Pro 4000 into the 9600 Plus Calibration Tester ear device port (large port). Place the probe firmly into the large port. See (Figure 54)“Braun Ear Device Port” on page 83. The temperature selection is non-specific for this particular test.
5. Wait approximately three seconds, press the thermometer **Start** button, and watch for the “Exac Temp” light to flash. See “Braun Exac Temp Light”.
6. Leave the thermometer in the 9600 Plus Calibration Tester until you hear a beep. Remove the Pro 4000 Thermometer from the 9600 Plus Calibration Tester and read the temperature in the thermometer's display.
7. Insert the Braun ThermoScan PRO 4000 thermometer into the UUT (should be on and in normal mode). The reading from the Braun ThermoScan PRO 4000 thermometer should transfer automatically to the UUT display. Verify this action.

## Changing the Calibration Set Point Temperature

To scroll from one set point to the next, press and hold the Temperature Selection button until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, starts to flash, and continues flashing until the cavity reaches the equilibrium at the new set point.

### To begin functional verification of the Masimo SpO2

Use the “Masimo Tester” to verify the Masimo SpO2 functionality.

**Figure 57. Masimo Tester**

1. With the UUT OFF, orient the Masimo Tester such that the mini-D connector mates with the SpO2 connector on the UUT; or use appropriate Masimo extension cable to connect the Masimo Tester and UUT.
2. Power ON the Spot LXi UUT and confirm the SpO2 reading in the Display Window is  $81\% \pm 3\%$  and the pulse reading is  $61 \text{ bpm} \pm 1 \text{ bpm}$ .

**Figure 58. SpO2 Reading**

3. Place the thumb and index finger on the gray buttons on either side of the Masimo Tester connector or extension cable connector; press the buttons firmly, and gently pull to remove the tester.

**Note** The functionality of the SpO2 module is not adjustable. If the SpO2 is not functioning properly, or the SpO2 and pulse performance cannot be verified against the calibrated simulator, contact Welch Allyn Technical Support [page iv](#).

### To begin functional verification of the Nellcor SpO2

Use a Nellcor-approved SpO2 simulator (SRC-MAX) to verify the Nellcor SpO2 functionality.

**Figure 59. Nellcor SpO2 Simulator**



1. With the UUT ON and completely booted, orient the simulator such that the connector mates with the SpO2 connector on the UUT; or use appropriate Nellcor extension cable to connect the simulator and UUT.

**Note** The UUT MUST be powered on prior to connecting the Nellcor SpO2 simulator.

2. On the Nellcor simulator change the SpO2 percentage from default 75% to 90%.
3. Confirm the SpO2 reading in the Display Window is  $90\% \pm 1\%$  and the pulse reading is  $60 \text{ bpm} \pm 1 \text{ bpm}$ .

**Figure 60. SpO2 Reading and Pulse Reading**

4. Remove the simulator by placing the thumb and index finger on either side of the simulator connector or extension cable connector; press the buttons firmly and gently pull to remove.

**Note** The functionality of the SpO2 module is not adjustable. If the SpO2 is not functioning properly, or the SpO2 and pulse performance cannot be verified against the calibrated simulator, contact Welch Allyn Technical Support [page iv](#)



# 7 Specifications

---

## Performance

This section describes normal ranges for Spot Vital Signs LXi.

### Blood Pressure Accuracy

Blood pressure accuracy meets or exceeds AAMI SP10:2002 standards for non-invasive blood pressure accuracy ( $\pm 5$  mmHg mean error, 8 mmHg standard deviation). Blood pressure accuracy is validated for pressure measurement using the upper arm only.

|                                                                     |                                                                                                                                                                                               |
|---------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Cuff Pressure Range</b>                                          | 0 to 300 mmHg                                                                                                                                                                                 |
| <b>Systolic Range</b>                                               | 60 to 250 mmHg                                                                                                                                                                                |
| <b>Diastolic Range</b>                                              | 30 to 160 mmHg                                                                                                                                                                                |
| <b>Blood Pressure Determination Time</b>                            | Typical: 15 seconds                                                                                                                                                                           |
| <b>Mean Arterial Pressure Range</b>                                 | 40 to 190 mmHg                                                                                                                                                                                |
| <b>Pulse Rate Range<br/>(using SpO<sub>2</sub> determination)</b>   | 25 to 240 bpm                                                                                                                                                                                 |
| <b>Pulse Rate Range<br/>(using Blood Pressure determination)</b>    | 35 to 199 bpm                                                                                                                                                                                 |
| <b>Pulse Rate Accuracy (using SpO<sub>2</sub> determination)</b>    | Without Motion: 25 to 245 bpm $\pm 3$ digits <sup>1</sup><br>With Motion: normal physiologic range (55 to 125 bpm) $\pm 5$ digits<br>Low Perfusion: 25 to 245 bpm $\pm 3$ digits <sup>1</sup> |
| <b>Pulse Rate Accuracy<br/>(using Blood Pressure determination)</b> | $\pm 5.0\%$                                                                                                                                                                                   |
| <b>Overpressure Cutoff</b>                                          | 315 mmHg $\pm 15$ mmHg                                                                                                                                                                        |

<sup>1</sup> Specification applies to device performance and was validated with Biotek and Nellcor simulators.

# Temperature Specifications

|                                                            |                                                                                                              |
|------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Temperature Range</b>                                   |                                                                                                              |
| SureTemp Plus                                              | 80° to 110° F (26.7° to 43.3° C)                                                                             |
| Braun ThermoScan PRO 4000                                  | 68° to 108° F (20° to 42.2° C)                                                                               |
| <b>Calibration Accuracy</b>                                |                                                                                                              |
| SureTemp Plus                                              | ±0.2° F (0.1° C) (Monitor Mode)                                                                              |
| Braun ThermoScan PRO 4000 for displayed temperature ranges | ± 0.4° F (±0.2° C) 95.9° to 107.6° F (35.5° to 42° C)<br>±0.5° F (±0.25° C) (outside this temperature range) |
| <b>Display Resolution</b>                                  | 0.1° F or ° C                                                                                                |
| <b>SureTemp Plus Predict Time</b>                          |                                                                                                              |
| Oral                                                       | Approx. 4 to 6 seconds                                                                                       |
| Adult Axillary                                             | Approx. 12 to 15 seconds (age 18 years and older)                                                            |
| Pediatric Axillary                                         | Approx. 10 to 13 seconds (age 17 years and younger)                                                          |
| Rectal                                                     | Approx. 10 to 13 seconds                                                                                     |

# SpO<sub>2</sub> Specifications

## Masimo Sensor Accuracy Guide

Accuracy specified when used with Masimo SET pulse oximetry monitors or with licensed Masimo SET pulse oximetry modules using PC series patient cables, during no motion. Numbers present  $\pm 1$  standard deviation. Plus or minus one standard deviation represents 68% of the population. SpO<sub>2</sub> accuracy from 70% to 100%. Pulse rate accuracy from 25 to 240 bpm.

### Performance Measurement Range

SpO<sub>2</sub>: 1 to 100%  
Pulse Rate: 25 - 240 beats per minute (BPM)

### Perfusion

0.02% to 20%

### Accuracy

Saturation: 70% to 100%  
No Motion: Adults, Pediatrics  $\pm 2$  digits  
Neonates  $\pm 3$  digits

Motion: Adults, Pediatrics  $\pm 3$  digits  
Neonates  $\pm 3$  digits

Low Perfusion: Adults, Pediatrics  $\pm 2$  digits  
Neonates  $\pm 3$  digits

### Pulse Rate Accuracy

Pulse Rate: 25 to 250 bpm

No Motion: Adults, Pediatrics, Neonates  $\pm 3$  digits

Motion: Adults, Pediatrics, Neonates  $\pm 5$  digits

Low Perfusion: Adults, Pediatrics, Neonates  $\pm 5$  digits

**Table 1. Masimo Sensor Accuracy Guide**

| Sensor     | Weight Range | Saturation Accuracy |           | Pulse Rate Accuracy |             |
|------------|--------------|---------------------|-----------|---------------------|-------------|
|            |              | No Motion           | Motion    | No Motion           | Motion      |
| LNOP-ADT   | > 30 kg      | $\pm 2\%$           | $\pm 3\%$ | $\pm 3$ bpm         | $\pm 5$ bpm |
| LNOP-PDT   | 10 to 50 kg  | $\pm 2\%$           | $\pm 3\%$ | $\pm 3$ bpm         | $\pm 5$ bpm |
| LNOP-NEO   | < 10 kg      | $\pm 3\%$           | $\pm 3\%$ | $\pm 3$ bpm         | $\pm 5$ bpm |
| LNOP-NEOPT | < 1 kg       | $\pm 3\%$           | $\pm 3\%$ | $\pm 3$ bpm         | $\pm 5$ bpm |
| LNOP-DCI   | > 30 kg      | $\pm 2\%$           | $\pm 3\%$ | $\pm 3$ bpm         | $\pm 5$ bpm |
| LNOP-DCIP  | 10 to 50 kg  | $\pm 2\%$           | $\pm 3\%$ | $\pm 3$ bpm         | $\pm 5$ bpm |

## Masimo Patents

The Masimo sensors and cables are covered under one or more of the following U.S.A. patents: 5,758,644; 5,823,950; 6,011,986; 6,157,850; 6,263,222; 6,501,975; and other applicable patents listed at [www.masimo.com/patents.htm](http://www.masimo.com/patents.htm).

## Nellcor Sensor Accuracy Guide

Accuracy specifications are based on controlled hypoxia studies with healthy, non-smoking adult volunteers over the specified saturation SpO<sub>2</sub> range. Pulse oximeter SpO<sub>2</sub> readings were compared to SaO<sub>2</sub> values of drawn blood samples measured by hemoximetry. All accuracies are expressed as  $\pm$  "X" digits. This variation equals  $\pm$  one standard deviation ( $\pm$  1 SD), which encompasses 68% of the population.

**Table 2. OxiMax Sensor Models, Single Patient Use**

| Sensor Models       | SpO <sub>2</sub> Range<br>70% to 100% |
|---------------------|---------------------------------------|
| MAX-AI <sup>1</sup> | $\pm$ 2                               |
| MAX-PI*             | $\pm$ 2                               |
| MAX-II              | $\pm$ 2                               |
| MAX-RI <sup>2</sup> | $\pm$ 3.5                             |

<sup>1</sup> The accuracy specification under motion conditions is  $\pm$  3. For a definition of motion, contact Nellcor Technical Services or your local Nellcor representative.

<sup>2</sup> The accuracy specification has been determined between saturations of 80% to 100%.

**Table 3. OxiCliq Sensor Models, Single Patient Use**

| Sensor Models | SpO <sub>2</sub> Range<br>70% to 100% |
|---------------|---------------------------------------|
| OXICLIQ-PI    | $\pm$ 2.5                             |

**Table 4. Reusable Sensor Models**

| Sensor Models              | SpO <sub>2</sub> Range<br>70% to 100% |
|----------------------------|---------------------------------------|
| D-YS (Infant to Adult)     | $\pm$ 3                               |
| D-YS and D-YSE             | $\pm$ 3.5                             |
| D-YS and D-YSPD            | $\pm$ 3.5                             |
| DS-100A                    | $\pm$ 3                               |
| OXI-A/N (Adult/neonate)    | Adult: $\pm$ 3<br>Neonate: $\pm$ 4    |
| OXI-P/I (Pediatric/infant) | $\pm$ 3                               |

## Nellcor Patents

Covered by one or more of the following U.S. patents and foreign equivalents:

5,485,847; 5,676,141; 5,743,263; 6,035,223; 6,226,539; 6,411,833; 6,463,310; 6,591,123; 6,708,049; 7,016,715; 7,039,538; 7,120,479; 7,120,480; 7,142,142; 7,162,288; 7,190,985; 7,194,293; 7,209,774; 7,212,847; 7,400,919.

## Mechanical

### Dimensions

Height: 10.63" (27 cm)  
 Length/Braun: 8" (20.32 cm)  
 Length/SureTemp Plus: 7.5" (19 cm)  
 Depth: 5.25" (13.34 cm)

### Weight

7.5 lbs (3.4 kg)

### Mounting

Self-supporting on rubber feet  
 Custom mobile stand  
 Custom wall mount

### Portability

May be hand-carried when held by the rear handle

## Electrical

### Power Requirements

Patient-rated isolation transformer is connected to AC mains:  
 100-240V, 50-60Hz, 0.4A or internal power: 6Vdc, 6Ah

## Patents

For patent information, please visit [www.welchallyn.com/patents](http://www.welchallyn.com/patents)

## Batteries

### Spot Vital Signs LXi

Sealed Lead Acid Battery (4500-84)  
 with external charger (4500-101A)

### Braun ThermoScan PRO 4000 thermometer

Rechargeable Battery Pack (53020-0000)

### External Printer

Ni-Mh Rechargeable Battery Pack (4500-505)  
 Includes a separate charger (713703 or BC1300M)

- The Spot Vital Signs LXi battery is 90 to 100% charged after 6 hours of charging. The Spot LXi battery automatically charges when Spot LXi is powered through the AC power transformer. An operator can use the device while the battery is charging; however, the battery charges faster when the instrument is not in operation.
- The rechargeable batteries in the Braun ThermoScan PRO 4000 thermometer require an additional 1 hour to charge.
- Charge the external printer battery 9 hours before initial use.

## Environmental



**WARNING** This device is not suitable for use in the presence of a flammable anesthetic mixture with air or oxygen or nitrous oxide. An explosion may result.

|                                      |                                                                                                                               |
|--------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Operating Temperature</b>         | 50° to 104° F (10° to 40° C)                                                                                                  |
| <b>Storage/Transport Temperature</b> | Device with SureTemp Plus: -13° to 131°F (-25° to 55°C)<br>Device with Braun ThermoScan PRO 4000: -4° to 122°F (-20° to 50°C) |
| <b>Relative Humidity</b>             | 15 to 95% (non-condensing)                                                                                                    |
| <b>Operating Altitude</b>            | -557 to 16,000 ft. (-170 to 4877 m)                                                                                           |

## Radio

## Mechanical

|                              |                                                                 |
|------------------------------|-----------------------------------------------------------------|
| <b>Dimensions</b>            | 3.93" L x 2.38" W x 1.42" T<br>99.78mm L x 60.5 mm W x 36.1mm T |
| <b>Mounting</b>              | Nylon Bracket                                                   |
| <b>Operating Temperature</b> | -40° C to +85° C                                                |
| <b>Storage Temperature</b>   | -40° C to +125° C                                               |

## Configuration Options

|                            |                                                                                                                                                                 |             |               |             |          |          |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|---------------|-------------|----------|----------|
| <b>DHCP Enabled or Not</b> | If not DHCP: <ul style="list-style-type: none"> <li>• IP Address</li> <li>• Subnet Address</li> <li>• Default Gateway Address</li> <li>• Primary DNS</li> </ul> |             |               |             |          |          |
| <b>Encryption</b>          | Wireless Security Mode                                                                                                                                          | Disabled    | WEP 64-bit    | WEP 128-bit | WPA2-PSK | WPA2-AES |
|                            | WPA2 Security (WPA2-PSK and WPA2-AES only)                                                                                                                      | CCMP        | CCMP and TKIP |             |          |          |
|                            | EAP Settings (WPA2-AES only)                                                                                                                                    | EAP Type    | TLS           | TTLS        | PEAP     | Auto     |
| <b>Channel</b>             |                                                                                                                                                                 |             |               |             |          |          |
| <b>Default WEP Key 1-4</b> |                                                                                                                                                                 |             |               |             |          |          |
| <b>SSID</b>                | WEP Authentication                                                                                                                                              |             |               |             |          |          |
| <b>WPA Pass Phrase</b>     |                                                                                                                                                                 |             |               |             |          |          |
| <b>WPA2 Pass Phrase</b>    | PSK Format                                                                                                                                                      | Hexadecimal | ASCII         |             |          |          |

# Guidance and Manufacturer's Declaration For Spot Vital Signs LXi Radio

## Emissions and Immunity Information

### Electromagnetic Emissions

The Spot Vital Signs LXi is intended for use in the electromagnetic environment specified below. The customer or user of the Spot Vital Signs LXi should assure that it is used in such an environment.

| <b>Emissions Test</b>                                       | <b>Compliance</b> | <b>Electromagnetic Environment - Guidance</b>                                                                                                                                                                                                                                                                                                                                                                                              |
|-------------------------------------------------------------|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions<br>CISPR 11                                    | Group 1           | The Spot Vital Signs LXi uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>The Spot Vital Signs LXi is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions<br>CISPR 11                                    | Class B           |                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Harmonic emissions<br>IEC 61000-3-2                         | Class A           |                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Voltage fluctuations/<br>flicker emissions<br>IEC 61000-3-3 | Complies          |                                                                                                                                                                                                                                                                                                                                                                                                                                            |

### Electromagnetic Immunity

The Spot Vital Signs LXi is intended for use in the electromagnetic environment specified below. The customer or user of the Spot Vital Signs LXi should assure that it is used in such an environment.

| <b>Immunity Test</b>                             | <b>IEC 60601 Test Level</b>                                  | <b>Compliance Level</b>                                      | <b>Electromagnetic Environment - Guidance</b>                                                                                                  |
|--------------------------------------------------|--------------------------------------------------------------|--------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ± 6 kV contact<br>± 8 kV air                                 | ± 6 kV contact<br>± 8 kV air                                 | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst<br>IEC 61000-4-4 | ±2 kV for power supply lines<br>±1 kV for input/output lines | ±2 kV for power supply lines<br>±1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment.                                                            |
| Surge<br>IEC 61000-4-5                           | ±1 kV differential mode<br>±2 kV common mode                 | ±1 kV differential mode<br>±2 kV common mode                 | Mains power quality should be that of a typical commercial or hospital environment.                                                            |

|                                                                                                              |                                                                                                            |                                                                                                            |                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Voltage dips, short interruptions, and voltage variations on power supply input lines.<br><br>IEC 61000-4-11 | >95% dip in 0.5 cycle<br><br>60% dip in 5 cycles<br><br>30% dip for 25 cycles<br><br>>95% dip in 5 seconds | >95% dip in 0.5 cycle<br><br>60% dip in 5 cycles<br><br>30% dip for 25 cycles<br><br>>95% dip in 5 seconds | Mains power quality should be that of a typical commercial or hospital environment. If the user of the Spot Vital Signs LXi requires continued operation during power mains interruptions, it is recommended that the Spot Vital Signs LXi be powered from an uninterruptible power supply or battery. |
| Power frequency (50/60Hz) magnetic field<br><br>IEC 61000-4-8                                                | 3 A/m                                                                                                      | 3 A/m                                                                                                      | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                              |

### Electromagnetic Immunity

The Spot Vital Signs LXi is intended for use in the electromagnetic environment specified below. The customer or user of the Spot Vital Signs LXi should assure that it is used in such an environment.

| Immunity Test                 | IEC 60601 Test Level        | Compliance Level | Electromagnetic Environment - Guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-------------------------------|-----------------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 kHz to 80 MHz | 3 Vrms           | Portable and mobile RF communications equipment should be used no closer to any part of the Spot Vital Signs LXi, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>Recommended separation distance</b><br>$d = (1.17) \sqrt{P}$                                                                                                                                                                                                                                                                                                                                        |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            | $d = (1.17) \sqrt{P}$ 80 MHz to 800 MHz<br>$d = (2.33) \sqrt{P}$ 800 MHz to 2.5 GHz<br><br>where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Spot Vital Signs LXi is used exceeds the applicable RF compliance level above, the electrocardiograph should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the electrocardiograph.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.



### Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Spot Vital Signs LXi

The Spot Vital Signs LXi is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Spot Vital Signs LXi can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Spot Vital Signs LXi as recommended below, according to the maximum output power of the communications equipment.

| Rated Max. Output Power of Transmitter (W) | Separation Distance According to Frequency of Transmitter (m) |                                            |                                             |
|--------------------------------------------|---------------------------------------------------------------|--------------------------------------------|---------------------------------------------|
|                                            | 150 kHz to 80 MHz<br>$d = (1.17) \sqrt{P}$                    | 80 MHz to 800 MHz<br>$d = (1.17) \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = (2.33) \sqrt{P}$ |
| 0.01                                       | 0.117                                                         | 0.117                                      | 0.233                                       |
| 0.1                                        | 0.37                                                          | 0.37                                       | 0.74                                        |
| 1                                          | 1.17                                                          | 1.17                                       | 2.33                                        |
| 10                                         | 3.70                                                          | 3.70                                       | 7.37                                        |
| 100                                        | 11.70                                                         | 11.70                                      | 23.30                                       |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.



## 8

# Troubleshooting

This chapter provides helpful information in troubleshooting the Spot Vital Signs LXi.

## Error Codes

The following tables of error codes provide a quick reference of the descriptions and probable causes of error codes.

**Table 1. Informational Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                                                                                                    | Corrective Action                                                              |
|--------------------------------|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|
| _n/a                           |                       | Copyright notice displayed on powers up splash screen.                                                                                            |                                                                                |
| _n/a                           |                       | Message displays when the user attempts to start another print while a print job is currently printing.                                           | Printing is already in progress                                                |
| _n/a                           |                       | Clear current reading.                                                                                                                            |                                                                                |
| _n/a                           |                       | Erase all readings.                                                                                                                               |                                                                                |
| _n/a                           |                       | Information screen telling the user how to operate the location identifier edit menu in Configuration Mode.                                       | Use the navigation keys to edit characters and the SELECT key to stop editing. |
| _n/a                           |                       | Message displays when Mod F cannot perform a STEP NIBP measurement and switches to normal method.                                                 | Restrict patient movement.                                                     |
| _n/a                           |                       | Message displays when in MEM mode and there are no readings to display.                                                                           |                                                                                |
| _n/a                           |                       | External devices configuration menu message when the user tries to exit the menu having enabling the information system, but not that Patient ID. | User must enable the Patient ID if using the information system.               |
| _n/a                           |                       | Message displays below txtPIDMsgB.                                                                                                                | User must press SELECT and correct the problem.                                |
| _n/a                           |                       | Information displays when the user enables device in external devices configuration menu.                                                         | User must connect the Spot Vital Signs LXi to an external device.              |
| _n/a                           |                       | Message displays below txtPIDMsgA                                                                                                                 | User must press SELECT when completed.                                         |
| _n/a                           |                       | Line displays below error/information message if the user can clear the error/condition.                                                          | User must press ACCEPT to clear.                                               |

**Table 1. Informational Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                                                          | Corrective Action                                                        |
|--------------------------------|-----------------------|---------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|
| _n/a                           |                       |                                                                                                         | User decides whether to reset everything back go factory default values. |
| _n/a                           |                       | Error condition occurs causing the unit to automatically shut down. The user cannot clear this message. | Spot Vital Signs LXi will shut down in 1 minute.                         |
| _n/a                           |                       |                                                                                                         | User confirms to start printing.                                         |

**Table 2. Blood Pressure Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description             | Corrective Action                                                                                             |
|--------------------------------|-----------------------|--------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| 10                             | C01                   | User abort.                                |                                                                                                               |
| 11                             | C03                   | Inflation too quick.                       | Check tubing and connections.                                                                                 |
| 12                             | E10                   | Overpressure                               | Check patient.                                                                                                |
| 13                             | C04                   | Major air leak.                            | Excessive inflation time; check for air leaks.                                                                |
| 14                             | C05                   | Excessive arm movement                     | Excessive noise, check cuff placement; restrict patient movement.                                             |
| 15                             | C07                   | Auto cycle skipped, system busy.           | Device will power down.                                                                                       |
| 16                             | C07                   | Auto cycle skipped, NIBP reached too soon. | Device will power down.                                                                                       |
| 17                             | C02                   | Auto zero failure.                         | Auto zero failure; check air obstruction; limit patient movement.                                             |
| 18                             | C04                   | Could not reach target pressure.           | Excessive inflation time; check for air leaks.                                                                |
| 19                             | E20                   | Missed A/D sample.                         | General internal malfunction. Call Technical Support.                                                         |
| 20                             | E20                   | A/D timeout failure.                       | General internal malfunction. Call Technical Support.                                                         |
| 21                             |                       | NIPB NOT AVAILABLE.                        | Correct probe tip conditions. If problem persists, replace probe. If problem still persists, replace module.  |
| 22                             | E11                   | Malfunction                                | User should try another cycle. If problem persists, replace probe.                                            |
| 23                             | C06                   | Reinflation needed.                        | Measurement was outside of device's measurement range.                                                        |
| 24                             | C05                   | Cycle took too long.                       | Check cuff placement; restrict patient movement.                                                              |
| 25                             | C12                   | Ambient temp out of range.                 | Check internal temperature at probe well. Verify ambient temperature does not exceed operating specification. |
| 26                             | C13                   | Battery failure                            | Battery failure; use wall transformer.                                                                        |

**Table 2. Blood Pressure Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                                                                                      | Corrective Action                                                                           |
|--------------------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
| 27                             | E11                   | > 10mm residual pressure > minutes.                                                                                                 | Verify that the supply voltage is within specification.                                     |
| 28                             | E11                   | >15mm residual pressure >3 mins.                                                                                                    | Verify that the supply voltage is within specification.                                     |
| 29                             | C02                   | Couldn't bleed pressure < 9 mmHg                                                                                                    | Auto zero failure, check tubing for kinks and connection integrity, limit patient movement. |
| 32                             |                       | Thermistor pulled away from the tip or heater broken.                                                                               | Try another predict cycle. If the problem persists, replace probe.                          |
| 33                             |                       | Probe not responsive. Probe not characterized/calibrated.                                                                           | Replace probe.                                                                              |
| 34                             |                       | Probe not characterized/calibrated.                                                                                                 | Replace probe.                                                                              |
| 41                             | C06                   | Low inflation; first step had largest value, first step's value > systolic value.                                                   | Measurement was outside of Spot Vital Signs LXi measurement range.                          |
| 42                             | C06                   | Dump too soon; last step had largest value, next to last step's value > diastolic value; will not normally occur in auto-dump mode. | Check connections; restrict patient movement.                                               |
| 43                             | C05                   | Step too big systolic; difference in pressure between step above and step below systolic is too great.                              | Excessive noise; check cuff placement and limit movement.                                   |
| 44                             | C05                   | Step too big diastolic; difference in pressure between step above and step below diastolic is too great.                            | Excessive noise; check cuff placement and limit movement.                                   |
| 45                             | C05                   | Range error: calculated BP is outside of measurement range.                                                                         | Excessive noise; check cuff placement and limit movement.                                   |
| 46                             | C05                   | Pressure reverse systolic; pressure of step above systolic is lower than pressure of step below systolic.                           | Excessive noise; check cuff placement and limit movement.                                   |
| 47                             | C05                   | Pressure reverse diastolic; pressure of step above systolic is lower than pressure of step below diastolic.                         | Excessive noise; check cuff placement and limit movement.                                   |
| 48                             | C05                   | Step too big MAP; difference in pressure between adjacent steps near MAP is too great.                                              | Excessive noise; check cuff placement and limit movement.                                   |
| 49                             | C05                   | Too few steps; less than 5 steps                                                                                                    | Excessive noise; check cuff placement and limit movement.                                   |
| 50                             | C05                   | Bad value diastolic; value of step below diastolic and step below that step are not consistent.                                     | Excessive noise; check cuff placement and limit movement.                                   |
| 51                             | C05                   | Value reverse systolic; value of step above systolic is greater than value of step below systolic.                                  | Excessive noise; check cuff placement and limit movement.                                   |
| 52                             | C05                   | Value reverse diastolic; value of step above systolic is greater than value of step below diastolic.                                | Excessive noise; check cuff placement and limit movement.                                   |

**Table 2. Blood Pressure Error Codes**

| <b>Internal Event Log Code Number</b> | <b>Displayed Device Code</b> | <b>Internal Event Log Description</b>                                                                                                                                                                                                                                      | <b>Corrective Action</b>                                  |
|---------------------------------------|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|
| 53                                    | C05                          | No info - can't be produced by Mod F if they don't exist in the code.                                                                                                                                                                                                      | Excessive noise; check cuff placement and limit movement. |
| 54                                    | C05                          | No info - can't be produced by Mod F if they don't exist in the code.                                                                                                                                                                                                      | Excessive noise; check cuff placement and limit movement. |
| 55                                    | C05                          | No info - can't be produced by Mod F if they don't exist in the code.                                                                                                                                                                                                      | Excessive noise; check cuff placement and limit movement. |
| 56                                    | C05                          | Out of memory; deflation contains too many steps. SW sets allowable limit of steps that can be stored and used for BP determination for a given cycle.                                                                                                                     | Excessive noise; check cuff placement and limit movement. |
| 57                                    | C05                          | Retarget error; cycle tried to retarget too many times. Only one allowed retarget allowed per cycle.                                                                                                                                                                       | Excessive noise; check cuff placement and limit movement. |
| 58                                    | C05                          | Repump error: Tried to repump too many times. Limit is set to 3 repumps. This can occur on a cycle with bad signal to noise ratio or if a subject's pulses are so small that pulses above the diastolic ratio are not detectable with the lowest allowed thread threshold. | Excessive noise; check cuff placement and limit movement. |
| 59                                    | C05                          | Systolic curve fit error; the steps on the systolic side of the envelope could not be curve fit. Causes include: too few or too many steps sent to curve fit routine; math error in curve fit routine.                                                                     | Excessive noise; check cuff placement and limit movement. |
| 60                                    | C05                          | Diastolic curve fit error; the steps on the systolic side of the envelope could not be curve fit. Causes include: too few or too many steps sent to curve fit routine; math error in curve fit routine.                                                                    | Excessive noise; check cuff placement and limit movement. |
| 61                                    | C05                          | Systolic slope error; slope of curve fit of systolic points has a slope out of range. Can occur with bad cuff fit, rapidly changing BP etc.                                                                                                                                | Excessive noise; check cuff placement and limit movement. |
| 62                                    | C05                          | Diastolic slope error; slope of curve fit of diastolic points has a slope out of range. Can occur with bad cuff fit, rapidly changing BP etc.                                                                                                                              | Excessive noise; check cuff placement and limit movement. |
| 63                                    | C05                          | Systolic extrapolation error; the systolic pressure that is represented by the systolic ratio is outside of the data sent to the curve fit routine.                                                                                                                        | Excessive noise; check cuff placement and limit movement. |
| 64                                    | C05                          | Diastolic extrapolation error; the systolic pressure that is represented by the diastolic ratio is outside of the data sent to the curve fit routine.                                                                                                                      | Excessive noise; check cuff placement and limit movement. |
| 65                                    | C05                          | Peak range error; the height of the peak of the envelope is outside normal limits.                                                                                                                                                                                         | Excessive noise; check cuff placement and limit movement. |
| 66                                    | C05                          | Reinflation too late.                                                                                                                                                                                                                                                      | Excessive noise; check cuff placement and limit movement. |

**Table 2. Blood Pressure Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                                          | Corrective Action                                                                                    |
|--------------------------------|-----------------------|-----------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
| 67                             | C06                   | High heart rate.                                                                        | Unable to determine blood pressure; pulse rate may be out-of-range.                                  |
| 68                             | Err                   | Too many reinflates.                                                                    | Restrict patient movement.                                                                           |
| 69                             | C06                   | Problem reading the module's EEPROM or saving to the module's EEPROM.                   | Verify that the supply voltage is correct. Try to make changes. If problem persists, replace module. |
| 70                             | C06                   | Problem reading the probe's EEPROM.                                                     | Verify that the supply voltage is correct. Try to make changes.                                      |
| 78                             | C06                   | The probe non-volatile error log memory error detection mechanism detected error.       | Measurement was outside of device's measurement range.                                               |
| 79                             | C06                   | The module's non-volatile error log memory error detecting mechanism detected an error. | Measurement was outside of device's measurement range.                                               |
| 5600                           | Err                   | Invalid ModF algorithm.                                                                 |                                                                                                      |
| 5601                           | Err                   | Invalid ModF init code.                                                                 |                                                                                                      |
| 5602                           | Err                   | Invalid ModF PMode.                                                                     |                                                                                                      |
| 5603                           | Err                   | Invalid ModF target pressure.                                                           |                                                                                                      |
| 5604                           | Err                   | Invalid ModF configuration.                                                             |                                                                                                      |

**Table 3. POST Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description | Corrective Action                                                                            |
|--------------------------------|-----------------------|--------------------------------|----------------------------------------------------------------------------------------------|
| 80                             | E30                   | FPROM checksum error.          | General internal malfunction; contact Technical Support.                                     |
| 81                             | E31                   | RAM test failure.              | General internal malfunction; contact Technical Support.                                     |
| 82                             | E32                   | Factory EEprom checksum error. | General internal malfunction; contact Technical Support.                                     |
| 83                             | E33                   | User EEprom checksum error.    | General internal malfunction; contact Technical Support.                                     |
| 84                             | E34                   | A/D converter failure.         | General internal malfunction; contact Technical Support.                                     |
| 85                             | E35                   | SpO2 module failure.           | General internal malfunction; contact Technical Support.                                     |
| 86                             | E36                   | Diatek module failure          | General internal malfunction; contact Technical Support.                                     |
| 87                             | E37                   | Printer problem.               | General internal malfunction; contact Technical Support.                                     |
| 88                             | E38                   | Real time clock failure.       | Date and time not set; see "Set the Date and Time Configurations" in the Directions for Use. |

**Table 3. POST Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description    | Corrective Action                                                                                                                                                               |
|--------------------------------|-----------------------|-----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 89                             | E39                   | No calibration signature.         | General internal malfunction; contact Technical Support.                                                                                                                        |
| 90                             | E40                   | Bad NIBP calibration record.      | General internal malfunction; contact Technical Support.                                                                                                                        |
| 91                             | E41                   | Stepper valve bad.                | General internal malfunction; contact Technical Support.                                                                                                                        |
| 92                             | E42                   | 8bit substitute for 16bit errors. | Internal communications error. Disconnect the battery and wait 5 minutes. Reconnect the battery, and then see "Set the Date and Time Configurations" in the Directions for Use. |
| 93                             | E43                   | Error log checksum error.         | General internal malfunction; contact Technical Support.                                                                                                                        |
| 94                             | Err-E44               | POST que full.                    | General internal malfunction; contact Technical Support.                                                                                                                        |
| 95                             | E44 E45               | POST que full.                    | General internal malfunction; contact Technical Support.                                                                                                                        |

**Table 4. SureTemp Plus Temperature Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                       | Corrective Action                                                 |
|--------------------------------|-----------------------|------------------------------------------------------|-------------------------------------------------------------------|
| 121                            | E0.4                  | Probe is over temperature (112° F / 44.4° C).        | Correct probe tip conditions. If problem persists, replace probe. |
| 122                            | E0.1                  | Excessive heater energy.                             | Try another predict cycle. If problem persists, replace probe.    |
| 123                            | 23                    | Host interface error.                                | Try another predict cycle. If problem persists, replace probe.    |
| 124                            | A^!                   | Device outside operating temperature range.          | Decrease ambient temperature.                                     |
| 125                            | Av!                   | Device outside operating temperature range.          | Increase ambient temperature.                                     |
| 126                            | 26                    | Invalid SureTemp algorithm.                          | Try another predict cycle. If problem persists, replace probe.    |
| 127                            | b^                    | Battery voltage over the algorithm's maximum value.  | Contact Technical Support.                                        |
| 128                            | bv                    | Battery voltage under the algorithm's minimum value. | Contact Technical Support.                                        |
| 129                            | 29                    | Battery voltage not set.                             | Try another predict cycle. If problem persists, replace probe.    |
| 130                            | 30                    | Prediction algorithm not set.                        | Try another predict cycle. If problem persists, replace probe.    |
| 131                            | 31                    | Ambient temperature not set.                         | Try another predict cycle. If problem persists, replace probe.    |



**Table 4. SureTemp Plus Temperature Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                           | Corrective Action                                                                                                                                                                                                                                     |
|--------------------------------|-----------------------|----------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 132                            | E0.5                  | No probe temperature rise.                               | Try another predict cycle. If problem persists, replace probe.                                                                                                                                                                                        |
| 133                            | E0.9                  | Probe EEPROM is uninitialized. Bad probe gain.           | Calibrate probe or replace probe.                                                                                                                                                                                                                     |
| 134                            | E0.9                  | Probe EEPROM is uninitialized. Bad probe response value. | Calibrate probe or replace probe.                                                                                                                                                                                                                     |
| 145                            | E0.2                  | A/D measurement overrun - pulse too long.                | Patient or environmental temperature conditions may be too low. Verify conditions at 80° F and 50° F respectively. If conditions are valid and problem persists, replace probe. If problem persists, replace the SureTemp Plus Thermometer module.    |
| 146                            | E0.2                  | A/D Measurement underrun - pulse too short.              | Patient or environmental temperature conditions may be too high. Verify conditions at 110° F and 104° F respectively. If conditions are valid and problem persists, replace probe. If problem persists, replace the SureTemp Plus Thermometer module. |
| 147                            | E4.1                  | A/D measurement overrun - RCAL pulse too long.           | If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                                                                                    |
| 148                            | E4.1                  | A/D measurement underrun - RCAL pulse too short.         | If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                                                                                    |
| 149                            | E4.0                  | A/D measurement overrun - PTB pulse too long.            | If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                                                                                    |
| 150                            | E4.0                  | A/D measurement underrun - PTB pulse too short.          | If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                                                                                    |
| 151                            | E0.2                  | A/D measurement overrun - timeout.                       | If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                                                                                    |
| 152                            | E0.9                  | Probe EEPROM is uninitialized.                           | Bad probe A/D calibration.                                                                                                                                                                                                                            |
| 159                            | bv                    | Displayed when battery voltage is too low.               |                                                                                                                                                                                                                                                       |
| 160                            | b^                    | Displayed when battery voltage is too high.              |                                                                                                                                                                                                                                                       |
| 161                            | E4.2                  | Reference power supply malfunction.                      | If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                                                                                    |
| 162                            | E4.8                  | Device not calibrated.                                   | Calibrate device.                                                                                                                                                                                                                                     |
| 163                            | E4.9                  | Missing probe well - cycle will not start.               | Insert probe well.                                                                                                                                                                                                                                    |
| 165                            | E4.3                  | Device EEPROM save failure.                              | Verify that the batteries are fresh. Try to make changes again. If problem persists, replace the SureTemp Plus Thermometer module.                                                                                                                    |
| 166                            | E0.8                  | EEPROM CRC Failure.<br>(CONFIG_TRANSFER_FAILURE)         | Verify that the batteries are fresh. Try again, if problem persists, replace the component that is indicated.                                                                                                                                         |

**Table 4. SureTemp Plus Temperature Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                        | Corrective Action                                                                                                                  |
|--------------------------------|-----------------------|-----------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| 167                            | E4.4                  | Device EEPROM read/save failure. (CONFIG_START_FAILURE)               | Verify that the batteries are fresh. Try to make changes again. If problem persists, replace the SureTemp Plus Thermometer module. |
| 168                            | E4.5                  | Device EEPROM read/save failure. (CONFIG_SLAVE_ACK_FAILURE)           | Verify that the batteries are fresh. Try to make changes again. If problem persists, replace the SureTemp Plus Thermometer module. |
| 169                            | E4.6                  | Device EEPROM read/save failure. (CONFIG_MASTER_ACK_FAILURE)          | Verify that the batteries are fresh. Try to make changes again. If problem persists, replace device.                               |
| 170                            | E0.6                  | Probe EEPROM read failure. (CONFIG_NULL_CHAR_FAILURE)                 | Verify that the batteries are fresh. Try to make changes again. If problem persists, replace the SureTemp Plus Thermometer module. |
| 171                            | 71                    | Device EEPROM access failure. (CONFIG_ACQUIRE_FAILURE)                | Try another predict cycle. If problem persists, replace probe.                                                                     |
| 172                            | 72                    | Device EEPROM access failure. (CONFIG_RELEASE_FAILURE)                | Try another predict cycle. If problem persists, replace probe.                                                                     |
| 173                            | 73                    | Device EEPROM access failure. (CONFIG_INVALID_PTR_FAILURE)            | Try another predict cycle. If problem persists, replace probe.                                                                     |
| 174                            | E4.7                  | Device EEPROM initialization failure. (CONFIG_INITIALIZE_FAILURE)     | If problem persists, replace the SureTemp Plus Thermometer module.                                                                 |
| 175                            | E0.7                  | Probe EEPROM is un-initialized. (CONFIG_UNINITIALIZED_PROBE)          | Verify that the batteries are fresh. If problem persists, replace probe.                                                           |
| 180                            | E5.0                  | Heater stuck on. (ERROR_HEATER_ON)                                    | If problem persists, replace the SureTemp Plus Thermometer module.                                                                 |
| 181                            | E5.0                  | Heater stuck off. (ERROR_HEATER_OFF)                                  | If problem persists, replace the SureTemp Plus Thermometer module.                                                                 |
| 182                            | E5.0                  | Heater on when should be off. (ERROR_HEATER_VOLTAGE_HIGH_Q)           | If problem persists, replace the SureTemp Plus Thermometer module.                                                                 |
| 183                            | E5.0                  | Heater on when should be off. (ERROR_HEATER_VOLTAGE_Q_HIGH_IMPEDANCE) | If problem persists, replace the SureTemp Plus Thermometer module.                                                                 |
| 184                            | E5.0                  | Heater is off when it should be on. (ERROR_HEATER_VOLTAGE_LOW)        | If problem persists, replace the SureTemp Plus Thermometer module.                                                                 |
| 185                            | E5.2                  | Heater hardware failsafe failure. (ERROR_HEATER_FAILSAFE)             | If problem persists, contact Technical Support.                                                                                    |

**Table 5. SpO<sub>2</sub> Error Codes**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description  | Corrective Action         |
|--------------------------------|-----------------------|---------------------------------|---------------------------|
| 1                              |                       | Missing or disconnected sensor. | Adjust or replace sensor. |
| 3                              | --                    | No pulse from sensor.           | Adjust or replace sensor. |

**Table 5. SpO<sub>2</sub> Error Codes**

| <b>Internal Event Log Code Number</b> | <b>Displayed Device Code</b> | <b>Internal Event Log Description</b>     | <b>Corrective Action</b>                   |
|---------------------------------------|------------------------------|-------------------------------------------|--------------------------------------------|
| 4                                     | E7                           | Communications failure.                   | Retake reading; contact Technical Support. |
| 5                                     |                              | 8-second timeout interval                 |                                            |
| 6                                     | C8                           | Defective or bad SpO <sub>2</sub> sensor. | Replace sensor.                            |
|                                       | C9                           | SpO <sub>2</sub> time limit exceeded.     | Remove sensor from patient.                |

**Table 6. Comms Errors**

| <b>Internal Event Log Code Number</b> | <b>Displayed Device Code</b> | <b>Internal Event Log Description</b>   | <b>Corrective Action</b>                   |
|---------------------------------------|------------------------------|-----------------------------------------|--------------------------------------------|
| 1000                                  | Err                          | Invalid Operation in this Context.      | Internal error, contact Technical Support. |
| 1001                                  | Err                          | Auto zero in progress.                  | Internal error, contact Technical Support. |
| 1002                                  | Err                          | Bad checksum.                           | Internal error, contact Technical Support. |
| 1003                                  | Err                          | Unimplemented command.                  | Internal error, contact Technical Support. |
| 1004                                  | Err                          | Buffer Overflow.                        | Internal error, contact Technical Support. |
| 1005                                  | Err                          | Overrun error.                          | Internal error, contact Technical Support. |
| 1006                                  | Err                          | Wrong data count.                       | Internal error, contact Technical Support. |
| 1007                                  | Err                          | FPROM erase error.                      | Internal error, contact Technical Support. |
| 1008                                  | Err                          | FPROM programming error                 | Internal error, contact Technical Support. |
| 1009                                  | Err                          | Bad Motorola record.                    | Internal error, contact Technical Support. |
| 1010                                  | Err                          | A/D channel locked.                     | Internal error, contact Technical Support. |
| 1011                                  | Err                          | RTC not present.                        | Internal error, contact Technical Support. |
| 1012                                  | Err                          | Attempt to reference nonexistent cycle. | Internal error, contact Technical Support. |
| 1013                                  | Err                          | Data range error.                       | Internal error, contact Technical Support. |
| 1014                                  | Err                          | Alarm pair range error.                 | Internal error, contact Technical Support. |
| 1015                                  | Err                          | No POST error to clear.                 | Internal error, contact Technical Support. |
| 1016                                  | Err                          | Can't clear this POST error.            | Internal error, contact Technical Support. |
| 1017                                  | Err                          | Command not command type.               | Internal error, contact Technical Support. |
| 1018                                  | Err                          | NAK from ST3800 protocol.               | Internal error, contact Technical Support. |
| 1019                                  | Err                          | No wireless link.                       | Internal error, contact Technical Support. |

**Table 7. Comms Host Errors**

| <b>Internal Event Log Code Number</b> | <b>Displayed Device Code</b> | <b>Internal Event Log Description</b> | <b>Corrective Action</b>                   |
|---------------------------------------|------------------------------|---------------------------------------|--------------------------------------------|
| 2001                                  | Err                          | Timeout.                              | Internal error, contact Technical Support. |
| 2002                                  | Err                          | Incorrect response header.            | Internal error, contact Technical Support. |
| 2003                                  | Err                          | Incorrect response checksum.          | Internal error, contact Technical Support. |
| 2004                                  | Err                          | Too much received data.               | Internal error, contact Technical Support. |
| 2005                                  | Err                          | User cancelled.                       | Internal error, contact Technical Support. |
| 2006                                  | Err                          | User paused.                          | Internal error, contact Technical Support. |

**Table 8. Printer Errors**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                         | Corrective Action                                                  |
|--------------------------------|-----------------------|------------------------------------------------------------------------|--------------------------------------------------------------------|
| 5400                           | Err                   | Problem with the printer                                               | Check all options; if problem persists, contact Technical Support. |
| 5401                           | Err                   | Printer is not responding.                                             |                                                                    |
| 5402                           | Err                   | Displayed when print started and printer is out of paper.              | Printer head will not shut all the way or too much paper.          |
| 5403                           | Err                   | Displayed when printing started and print head is not closed in place. | Add paper.                                                         |

**Table 9. Spot Vital Signs LXi Errors**

| Internal Event Log Code Number | Displayed Device Code | Internal Event Log Description                                    | Corrective Action                     |
|--------------------------------|-----------------------|-------------------------------------------------------------------|---------------------------------------|
| 5500                           |                       | Unknown error code.                                               |                                       |
| 5501                           | Err                   | Displayed when battery is too low for unit to continue operating. | Charge battery; use wall transformer. |
| 5502                           | Err                   | Key stuck error causing the unit to be automatically shut down.   |                                       |

**Table 10. 802.11 b Radio Errors**

| LED   | LED Color      | Function                                                                                                                                                        |
|-------|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power | Off            | Bridge is not receiving power.                                                                                                                                  |
|       | Red            | Bridge failed its Power On Self Test (POST) and is not configured for wireless communication.                                                                   |
|       | Amber          | Bridge passed its POST but is not configured for wireless communication.                                                                                        |
|       | Green          | Bridge passed its POST and is configured for wireless communication.                                                                                            |
| Link  | Off            | Bridge is not receiving power.                                                                                                                                  |
|       | Red            | Bridge is searching for an Access Point.                                                                                                                        |
|       | Green          | Wireless network and MAC have associated with an Access Point.                                                                                                  |
| Comm  | Off            | No power or no wireless TCP session is established and no physical serial connection is detected.                                                               |
|       | Red            | No wireless TCP session is established; a physical connection is detected.                                                                                      |
|       | Blinking Red   | A physical serial connection was detected and there is serial traffic present on that connection, but no wireless TCP session is established.                   |
|       | Amber          | A wireless TCP connection is established but no physical serial connection is detected (i.e. no serial cable is attached to the Bridge).                        |
|       | Blinking Green | A wireless TCP session is established, a physical serial connection is detected, and the Bridge is transmitting or receiving data across the wired serial port. |
|       | Green          | A wireless TCP session is established, a physical serial connection is detected, but there is no active data movement across the wired serial port.             |

**NOTE:**The 802.11a/b/g radio has no functional LED displays.

## Radio Test

Refer to Manufacturer's Directions for Use.

# Service Work Checklist

## Spot Vital Signs

| Model Number | Serial Number | BP Cycle Count | Technician | Date |
|--------------|---------------|----------------|------------|------|
|              |               |                |            |      |

| Test                                   | Test Data    | Pass/Fail | Test Specification             |
|----------------------------------------|--------------|-----------|--------------------------------|
| Unit SW Version                        |              | N/A       | N/A                            |
| SP02 SW Version                        |              | N/A       | N/A                            |
| Thermometry SW Version                 |              | N/A       | N/A                            |
| Unit Pressure @ 0 mmHg                 |              |           | 0 mmHg (+/- 1 mmHg)            |
| Unit Pressure @ 50 mmHg                |              |           | 50 mmHg (+/- 1 mmHg)           |
| Unit Pressure @ 150 mmHg               |              |           | 150 mmHg (+/- 1.5 mmHg)        |
| Unit Pressure @ 250 mmHg               |              |           | +/- 2.0 mmHg                   |
| Test Voltage @ 5.5 V                   |              |           | +0.1 VDC                       |
| Current Test – Sleep State             |              |           | < 0.5 A                        |
| Current Test – Patient DAQ             |              |           | < 1.25 A                       |
| Current Test – Valve Pump              |              |           | < 2.0 A                        |
| Noise Level                            |              |           | < 0.05 mmHg                    |
| Button Test                            |              |           | Pass/Fail                      |
| Interface Test                         |              |           | Pass/Fail                      |
| Leak Test                              |              |           | <= 6mmHg in 15 seconds         |
| Dump Test                              |              |           | < 10 Seconds                   |
| Inflation Test                         |              |           | < 7 Seconds                    |
| Valve Control Test                     |              |           | Pass/Fail                      |
| Over Pressure Test                     |              |           | 296 to 329 mmHg                |
| SP02 Sensor - Masimo                   |              |           | 81% +/- 3% 61bpm +/- 1bpm      |
| SP02 Sensor - Nellcor                  |              |           | 90% +/- 1% 60bpm +/- 1bpm      |
| Temperature Calibration Suretemp Plus  |              |           | 96.8 +/- 0.2°F 105.8 +/- 0.2°F |
| Temperature Calibration Braun Pro 4000 |              |           | 96.8 +/- 0.2°F 105.8 +/- 0.2°F |
| Temperature Cal Key Suretemp Plus      |              |           | 97.3 +/- 0.2°F                 |
| Radio Ping Test                        | Pass or Fail |           | Pass                           |

## Spot LXi Radio Configuration Tool

This CD configures the Spot Vital Signs LXi Radio to a local network. See the Welch Allyn Spot LXi Radio Configuration Procedure (80012309) for further information.



## 9

## Supplies and Accessories

## Blood Pressure

**Table 1. Reusable Blood Pressure Cuffs (1 per pack)**

| REF     | Size  | REF     | Size        |
|---------|-------|---------|-------------|
| 4500-01 | Child | 4500-03 | Large Adult |
| 4500-02 | Adult | 4500-04 | Thigh       |

**Table 2. Reusable Blood Pressure Cuffs (1 per pack)**

| REF          | Size        | REF          | Size        |
|--------------|-------------|--------------|-------------|
| REUSE-07-2MQ | Infant      | REUSE-11-2MQ | Adult       |
| REUSE-08-2MQ | Small Child | REUSE-12-2MQ | Large Adult |
| REUSE-09-2MQ | Child       | REUSE-13-2MQ | Thigh       |
| REUSE-10-2MQ | Small Adult |              |             |

**Table 3. Disposable Blood Pressure Cuffs (5 per box)**

| REF         | Size        | REF         | Size        |
|-------------|-------------|-------------|-------------|
| SOFT-07-2MQ | Infant      | SOFT-11-2MQ | Adult       |
| SOFT-08-2MQ | Small Child | SOFT-12-2MQ | Large Adult |
| SOFT-09-2MQ | Child       | SOFT-13-2MQ | Thigh       |
| SOFT-10-2MQ | Small Adult |             |             |

**Table 4. Miscellaneous Supplies and Accessories**

| REF     | Description                    | REF     | Description                   |
|---------|--------------------------------|---------|-------------------------------|
| 4500-30 | Blood Pressure Hose (5ft/1.5M) | 4500-31 | Blood Pressure Hose (10ft/3M) |
| 4500-32 | Blood Pressure Hose (8ft/2.4M) |         |                               |

# Temperature

**Table 5. SureTemp Plus**

| REF        | Description                                                           |
|------------|-----------------------------------------------------------------------|
| 02895-000  | SureTemp Plus Oral Probe and Well (9 feet/2.7M)                       |
| 02895-100  | SureTemp Plus Rectal Probe and Well (9 feet/2.7M)                     |
| 02894-0000 | SureTemp Plus Oral Well                                               |
| 02894-1000 | SureTemp Plus Rectal Well                                             |
| 06138-000  | SureTemp Plus Temperature Calibration Key                             |
| 01802-110  | 9600 Plus Calibration Tester                                          |
| 05031-101  | Disposable SureTemp Plus Probe Covers (1,000 covers, packaged 25/box) |
| 20500-025M | Probe Cover-Box 25                                                    |
| 400194     | Assy, Suretemp Plus Components, Spot Vital Signs LXi                  |
| 400195     | PCA, Suretemp Plus OEM without EDAC Calibrated                        |
| 400196     | PCA, EDAC, Spot Vital signs LXi                                       |

**Table 6. Braun ThermoScan PRO 4000**

| REF        | Description                                                                                                               |
|------------|---------------------------------------------------------------------------------------------------------------------------|
| 04000-200  | Braun ThermoScan PRO 4000 Thermometer (for North America, South America, and Asia Pacific)                                |
| 04000-600  | Braun ThermoScan PRO 4000 Thermometer (for Europe, Middle East, and Africa)                                               |
| 05075-800  | Braun ThermoScan PRO 4000 Disposable Probe Covers (Case of 800 covers for North America, South America, and Asia Pacific) |
| 04000-800  | Braun ThermoScan PRO 4000 Disposable Probe Covers (Case of 800 covers for Europe, Middle East, and Africa)                |
| 01802-110  | 9600 Plus Calibration Tester                                                                                              |
| 53020-0000 | Braun ThermoScan PRO 4000 Rechargeable Battery Pack                                                                       |
| 4500-53    | Braun Unlocking Key                                                                                                       |
| 4500-54    | Repair-Braun with rechargeable batteries                                                                                  |

# Pulse Oximetry

## Masimo Accessories

**Table 7. Adhesive Sensors: Single-Patient Use**

| Catalog #    | Description                                      | Weight Range                   | Quantity   |
|--------------|--------------------------------------------------|--------------------------------|------------|
| NEO-WRAP-RP  | Replacement Wrap Neo Adhesive                    |                                |            |
| LNCS-INF-3   | Disposable Adhesive Finger Sensor-Infant         |                                |            |
| INF-WRAP-RP  | Replacement Wrap Infant Adhesive                 |                                |            |
| YI-AD        | Multisite Adhesive Wrap AD/PED/NEO, YI Sensor    |                                |            |
| YI-FM        | Multisite Foam Wrap AD/PED/NEO, YI Sensor        |                                |            |
| LNCS-YI      | Multisite Reusable Sensor                        |                                |            |
| LNCS-TC-I    | Reusable Ear Sensor                              |                                |            |
| LNCS-NEO-L-3 | Disposable Adhesive Finger Sensor- NEO/Adult     |                                |            |
| LNCS-ADTX    | Adhesive Finger Sensor - Adult (20 per case)     | >30 kg                         | LNCS-ADTX  |
| LNCS-PDTX    | Adhesive Finger Sensor - Pediatric (20 per case) | 10 to 50 kg                    | LNCS-PDTX  |
| LNCS INF-L   | Adhesive Finger Sensor - Infant (20 per case)    | 3 to 20 kg                     | LNCS INF-L |
| LNOP-ADT     | Adult sensor                                     | >66 lbs (30 kg)                | 20         |
| LNOP-PDT     | Pediatric sensor                                 | 22 to 110 lbs<br>(10 to 50 kg) | 20         |
| LNOP-NEO     | Neonatal sensor                                  | <22 lbs (10 kg)                | 20         |
| LNOP-NEOPT   | SoftTouch neonatal preterm sensor                | <2.2 lbs (1 kg)                | 20         |

**Table 8. Reusable Sensor**

| Catalog # | Description                   | Weight Range    | Quantity |
|-----------|-------------------------------|-----------------|----------|
| LNCS-TC-I | Reusable Ear Sensor           |                 |          |
| LNCS-YI   | Multisite Reusable Sensor     |                 |          |
| LNCS-DCI  | Finger sensor - adult         | >66 lbs (30 kg) | 1        |
| LNCS-DCIP | Finger sensor - pediatric     | 10 to 50 kg     | 1        |
| LNOP-DCI  | Finger clip probe - adult     | >66 lbs (30 kg) | 1        |
| LNOP-DCIP | Finger clip probe - pediatric | 10 to 50 kg     | 1        |

**Table 9. Sensor Cables**

| Catalog # | Description                                | Weight Range | Quantity |
|-----------|--------------------------------------------|--------------|----------|
| LNC-4-WA  | 4-foot cable with DB-9 connector for LNCS  | NA           | 1        |
| LNC-10-WA | 10-foot cable with DB-9 connector for LNCS | NA           | 1        |
| PC-04     | 4-foot cable with sensor connector         | NA           | 1        |

**Table 9. Sensor Cables**

| Catalog # | Description                         | Weight Range | Quantity |
|-----------|-------------------------------------|--------------|----------|
| PC-08     | 8-foot cable with sensor connector  | NA           | 1        |
| PC-12     | 12-foot cable with sensor connector | NA           | 1        |

## Nellcor Accessories

**Table 10. OxiMax Adhesive Sensors: Single-Patient Use**

| Catalog # | Description                        | Weight Range                   | Quantity   |
|-----------|------------------------------------|--------------------------------|------------|
| MAX-AI    | MAX-A Adhesive Sensor, adult       | >66 lbs (30 kg)                | Case of 24 |
| MAX-PI    | MAX-P Adhesive Sensor, pediatric   | 22 to 110 lbs<br>(10 to 50 kg) | Case of 24 |
| MAX-II    | MAX-I Adhesive Sensor, infant      | 6.5 to 44 lbs<br>(3 to 20 kg)  | Case of 24 |
| MAX-RI    | MAX-R Adhesive Sensor, adult nasal | >110 lbs (50 kg)               | Case of 24 |

**Table 11. OxiMax OxiCliq® Sensors: Reusable Cable with Adhesive Sensor Bandage**

| Catalog # | Description                        | Weight Range                   | Quantity   |
|-----------|------------------------------------|--------------------------------|------------|
| OC-3      | OxiCliq Sensor Cable (3 ft / 91cm) | N/A                            | 1          |
| OXCIIQ-PI | OxiCliq P, pediatric               | 22 to 110 lbs<br>(10 to 50 kg) | Case of 24 |

**Table 12. OxiMax Reusable Sensors**

| Catalog # | Description                                   | Weight Range                               | Quantity |
|-----------|-----------------------------------------------|--------------------------------------------|----------|
| DS-100A   | Durasensor® DS-100A finger-clip sensor, adult | >88 lbs (40 kg)                            | 1        |
| OXI-A/N   | Oxiband® OXI-A/N, adult/neonatal              | < 6.5 lbs or > 88 lbs<br>(<3 kg or >40 kg) | 1        |
| OXI-P/I   | Oxiband OXI-P/I, pediatric/infant             | 6.5 lbs to 88 lbs<br>(3 to 40 kg)          | 1        |
| D-YS      | Dura-Y® D-YS, multisite sensor                | >2.2 lbs (1 kg)                            | 1        |
| D-YSE     | D-YSE ear clip for Dura-Y sensor              | >66 lbs (30 kg)                            | 1        |
| D-YSPD    | PediCheck™ D-YSPD pediatric spot-check sensor | 6.5lbs to 88 lbs<br>(3 to 40 kg)           | 1        |

**Table 13. OxiMax Sensor Cables**

| Catalog # | Description       | Weight Range | Quantity |
|-----------|-------------------|--------------|----------|
| DOC-10    | DOC-10 (10 ft/3M) | N/A          | 1        |

## Miscellaneous

| REF       | Description                                                        |
|-----------|--------------------------------------------------------------------|
| 4500-60   | Mobile Stand, Spot Vital Signs LXi                                 |
| 4500-62   | Wall Mount, Spot Vital Signs LXi                                   |
| 4500-84   | Lead Acid Battery                                                  |
| 705310    | Directions for Use, Spot Vital Signs LXi English                   |
| 705298    | Quick Reference Card                                               |
| 4500-11   | Child Inf. Bag, 2 Tube Monitor                                     |
| 4500-12   | Adult Bag, 2 Tube Monitor                                          |
| 4500-13   | Large Adult Inf. Bag, 2 Tube Monitor                               |
| 4500-14   | Thigh Inf. Bag, 2 Tube Monitor                                     |
| 4500-75   | Kit, Desk Plate, Spot Vital Signs LXi                              |
| 4500-89   | Service Manual, Spot Vital Signs LXi                               |
| 4500-150E | In-Service CD, Spot Vital Signs LXi                                |
| 4500-100  | Carrying Case, Spot Vital Signs LXi                                |
| 4500-101A | AC Power Transformer (desktop transformer, line cord not included) |
| 4500-200  | Cuffkit, Child, Adult, Thigh, 2T Cuff & Bag                        |
| 4500-202  | Cuffkit, Small Child, Child, Small Adult 2T                        |
| 4500-400  | Line Cord (United States/Canadian/Japanese version)                |
| 4500-402  | Line Cord (European version)                                       |
| 4500-404  | Line Cord (United Kingdom version)                                 |
| 4500-406  | Line Cord (Australian version)                                     |
| 4500-408  | Line Cord (South African version)                                  |
| 4500-450  | Transformer Mount, Spot Vital Signs LXi                            |
| 4500-500  | Kit, Printer, Spot Vital Signs LXi                                 |
| 4500-505  | Printer Rechargeable Replacement Battery                           |
| 4500-510  | Printer Paper (10 rolls)                                           |
| 4500-520  | Label Paper (10 rolls) and Cleaning Kit                            |
| 4500-900  | Kit, Dev Conn SDK Developers                                       |
| 4500-905  | Repair Tool Software CD Kit                                        |
| 4500-906  | Spot Vital Signs LXi Configuration CD                              |
| 4500-907  | Spot Vital Signs LXi Firmware Loader CD                            |
| 4500-915  | Image Team 4600 2D Bar Code Scanner Kit                            |
| 4500-920  | 802.11b (DEPAC) Radio Kit                                          |
| 4500-921  | Radio Directions for Use                                           |
| 4500-922  | 802.11 a/b/g Spot Vital Signs LXi Radio Kit                        |
| 4500-923  | Spot Vital Signs LXi a/b/g Radio Directions for Use                |

| REF      | Description                                                   |
|----------|---------------------------------------------------------------|
| 4500-925 | USB 2.0 Cable for Wired Connectivity                          |
| 4500-926 | Cable for Wired Connectivity, Keyspan (USB to serial adapter) |
| 4500-927 | Spot Vital Signs LXi USB/serial Cable Kit                     |
| 5200-08  | Calibration T Fitting                                         |

## Service Contracts

See the Welch Allyn Partners In Care™ web site (<http://www.welchallyn.com/promotions/services/default.htm>) for a description of the different service contracts offered.





## 10

## Field Replaceable Units

The following list identifies the available FRUs for Spot Vital Signs LXi. To order an FRU, contact Welch Allyn Technical Support.

**Table 1. Spot Vital Signs LXi Base Components**

| Component Number | Object Description                                  |
|------------------|-----------------------------------------------------|
| 400163           | ASSY, POWER HARNESS, Spot Vital Signs LXi           |
| 400170           | ASSY, PUMP HARNESS                                  |
| 403912           | PCA, MAIN, Spot Vital Signs LXi                     |
| 400200           | BASE UNIT, SURETEMP, NELLCOR - Spot Vital Signs LXi |
| 400203           | BASE UNIT COMPONENTS, Spot Vital Signs LXi          |
| 400207           | LANGUAGE COMPONENTS, ENG, Spot Vital Signs LXi      |
| 400387           | PBCA BP OEM, (MOD F)                                |
| 4500-901         | Spot Vital Signs LXi Pump and Pump Tubing Kit       |
| 4500-902         | Spot Vital Signs LXi Pump Tubing Kit                |
| 4500-84          | Lead Acid Battery                                   |
| 4500-101A        | AC Power Transformer                                |
| 400720           | ASSY, BALLAST HARNESS, Spot Vital Signs LXi         |
| 400732           | ASSY, BATTERY, Spot Vital Signs LXi                 |
| 401082           | Assy, Bezel and Window, Spot Vital Signs LXi        |
| 401250           | Assy, Fan Unit Cooling                              |
| 421051-12        | TUBING, 1/8 X 1/4 X 2.50                            |
| 421051-9         | TUBING, 1/8 X 1/4 X 0.95                            |
| 59P586           | FOOT, RUBBER                                        |
| 620028-E         | LABEL, CAUTION                                      |
| 620217           | FITTING, "T", 1/8 X 1/16X 1/8                       |
| 700101           | HSG, HANDLE, LXi                                    |
| 700102           | HSG, HANDLE INSERT, LXi Ultra                       |
| 700103           | HSG, REAR, LXi                                      |
| 700105           | HSG, BATT DOOR, Spot Vital Signs LXi                |
| 700110           | BUMPER, LCD, Spot Vital Signs LXi                   |
| 707237           | SWITCH ARRAY, Spot Vital Signs LXi-New              |
| 700119           | INVERTER, LCD                                       |

**Table 1. Spot Vital Signs LXi Base Components**

|            |                                                |
|------------|------------------------------------------------|
| 705653     | HOUSING, FRONT, Spot Vital Signs LXi           |
| 704423     | HSG, PANEL, SPO2-BLANK, Spot Vital Signs LXi   |
| 4500-89    | SERVICE MANUAL, Spot Vital Signs LXi           |
| 401427     | LABEL, ICON ID, Spot Vital Signs LXi           |
| 700863     | CBL,FLX, 1/2mm,15POS,1.0"                      |
| 704375     | FLEX, MAIN, Spot Vital Signs LXi               |
| 701040     | FITTING, PLASTIC L, 1/16 X 1/16                |
| 701502     | Valve,Pneumatic,6PSI,5V,TH                     |
| 703077     | STANDOFF, PCB, BROACHING, #4 X .500 LG         |
| 703118     | BATTERY CLIP, AA, SINGLE                       |
| 703354     | DISPLAY, LCD, 320X240, BLITE, MONO, 3.3V, 5.7" |
| 703355     | FILTER, AIR, 1/8 BARBS                         |
| 705284     | Spot Vital Signs LXi Warranty Card English     |
| 705298     | QRC, Spot Vital Signs LXi English              |
| 705310     | DFU, Spot Vital Signs LXi English              |
| 703842     | TUBING, SILICONE, .063ID X .125OD X 5.50       |
| 700847     | TUBING, SILICONE, .125IDx.250Dx1.75LG          |
| 421023-501 | ASSEMBLY, SPOT PUMP REPAIR                     |
| 4500-500   | KIT, PRINTER, SPOT ULTRA (LXi)                 |
| 4500-501   | Spot LXi PRINTER UPGRADE KIT                   |
| 713703     | SPOT LXi PRINTER POWER SUPPLY                  |

**Table 2. SureTemp Plus Components**

| Component Number | Object Description                                |
|------------------|---------------------------------------------------|
| 02692-100        | PROBE ASSY,ORAL 9FT                               |
| 02894-0000       | ASSY,PRB WLL,BLUE,OEM,M690_692                    |
| 400194           | ASSY, SURETEMP+, COMPONENTS, Spot Vital Signs LXi |
| 400195           | PCA, SURETEMP+, OEM W/O EDAC, CALIBRATED          |
| 400196           | PCA, EDAC, Spot Vital Signs LXi                   |
| 700111           | HOUSING, THERM OEM POD, Spot Vital Signs LXi      |
| 700112           | PLATE, POD MOUNTING, Spot Vital Signs LXi         |
| 700858           | CBL,FLX,1MM,6POS                                  |
| 700859           | CBL,FLX,1MM,8POS                                  |
| 70999-0000       | LABEL,SURETEMP TECHNOLOGY                         |
| 02692-100        | PROBE ASSY,ORAL 9FT                               |

**Table 3. Braun ThermoScan PRO 4000 Components**

| <b>Component Number</b> | <b>Object Description</b>                      |
|-------------------------|------------------------------------------------|
| 400652                  | PCA, BRAUN INTERFACE, Spot Vital Signs LXi     |
| 400831                  | BRAUN COMPONENTS, Spot Vital Signs LXi         |
| 700112                  | PLATE, POD MOUNTING, Spot Vital Signs LXi      |
| 700820                  | HSG, BOTTOM, BRAUN POD, Spot Vital Signs LXi   |
| 700822                  | HSG, BRAUN POD, Spot Vital Signs LXi           |
| 704365                  | WINDOW, BRAUN POD, Spot Vital Signs LXi        |
| 700859                  | CBL,FLX,1MM,8POS                               |
| 703161                  | SECURITY LOCK, BRAUN POD, Spot Vital Signs LXi |
| 703162                  | LATCH, BRAUN POD, Spot Vital Signs LXi         |
| 703164                  | CONTACT, BRAUN POGO, Spot Vital Signs LXi      |
| 703165                  | SPRING, COMPRESSION                            |
| 703192                  | COVER, BRAUN SECURITY, Spot Vital Signs LXi    |

**Table 4. Masimo SpO<sub>2</sub> Components**

| <b>Component Number</b> | <b>Object Description</b>                    |
|-------------------------|----------------------------------------------|
| 008-0648-00             | MASIMO LNCP-DCI ADULT REUSABLE SENSOR        |
| 008-0692-01             | MASIMO - 8' INTERFACE CABLE                  |
| 008-0960-00             | MASIMO-LNCS-DCI ADULT REUSABLE SENSOR        |
| 008-1013-00             | SPO2 CABLE, LNCS, MASIMO, MINI-D, 4', EXP    |
| 400192                  | PCA, MASIMO OEM, Spot Vital Signs LXi        |
| 400205                  | MASIMO SPO2 COMPONENTS, Spot Vital Signs LXi |
| 400555                  | PCA, MASIMO ADAPTOR, Spot Vital Signs LXi    |
| 400610                  | ASSY, CABLE, FLAT FLEX, SPO2                 |
| 704422                  | HSG, PANEL, SPO2-MASIMO, LXi                 |
| 705651                  | HOUSING, FRONT MASIMO, Spot LXi              |
| 711350                  | MASIMO MS-11 SPO2, STD, SW V4.6.0.2          |
| LNC-10-WA               | MASIMO, 10' EXT CABLE, LNCS/DB-9             |
| LNC-4-WA                | MASIMO, 4' EXT CABLE, LNC/DB-9               |
| LNCS-ADTX               | MASIMO, LNCS-ADTX ADULT ADHSV, 20 CASE       |
| LNCS-DCI                | MASIMO, LNCS-DCI ADULT REUSABLE SENSOR       |
| LNCS-DCIP               | MASIMO, LNCS-DCIP PED REUSABLE SENSOR        |
| LNCS-INF-L              | MASIMO, LNCS-INF-L INFANT ADHSV, 20 CASE     |
| LNCS-PDTX               | MASIMO, LNCS-PDTX PED ADHSV, 20 CASE         |
| LNOP-ADT                | DISPOSABLE ADHESIVE FINGER CLIP PROBE-AD     |
| LNOP-DCI                | REUSABLE FINGER CLIP PROBE-ADULT             |
| LNOP-DCIP               | REUSABLE FINGER CLIP PROBE-PEDIATRIC         |
| LNOP-NEO                | DISP ADHESIVE FINGER CLIP PROBE-NEONATE      |

**Table 4. Masimo SpO<sub>2</sub> Components**

| Component Number | Object Description                        |
|------------------|-------------------------------------------|
| LNOP-NEOPT       | DISPOSABLE FINGER SENSOR-NEONATE          |
| LNOP-PDT         | DISP ADHESIVE FINGER CLIP PROBE-PEDIATRIC |
| PC-04            | 4-FOOT CABLE W/ SENSOR CONN               |
| PC-08            | 8-FOOT CABLE W/ SENSOR CONN               |
| PC-12            | 12-FOOT CABLE W/ SENSOR CONN              |

**Table 5. Nellcor SpO<sub>2</sub> Components**

| Component Number | Object Description                           |
|------------------|----------------------------------------------|
| 400183           | PCA, NELCOR OEM, Spot Vital Signs LXi        |
| 400201           | NELCOR SPO2 COMPONENTS, Spot Vital Signs LXi |
| 400610           | ASSY, CABLE, FLAT FLEX, SPO2                 |
| 620377-1         | "NELLCOR WORKS HERE" LABEL                   |
| 704421           | HSG, PANEL, SPO2-NELLCOR, LXi                |
| 704870           | EVSPNELLCOR NELL3 PULSE OXIMETRY             |
| DOC-10           | CABLE, SPO2 EXTENSION, NELLCOR               |
| DEC-4            | EXTENSION CABLE 4 FT.                        |
| DS-100A          | DURASENSOR, ADULT, DS-100A                   |
| D-YS             | DURA-Y OXYGEN TRANSDUCER (1 sensor-40 wraps) |
| D-YSE            | EAR CLIP (USE WITH DURA-Y SENSOR)            |
| D-YSPD           | PEDICHECK PEDIATRIC SPOT CHECK               |
| MAX-AI           | NELLCOR OXIMAX ADULT SENSOR (SU CS/24)       |
| MAX-II           | NELLCOR OXIMAX INFANT SENSOR (SU CS/24)      |
| MAX-PI           | NELLCOR OXIMAX PED SENSOR (SU CS/24)         |
| MAX-RI           | NELLCOR OXIMAX ADULT NASAL SEN (CS/24)       |
| OC-3             | NELLCOR OXICLIQ SENSOR CABLE                 |
| OXI-A/N          | NELLCOR OXIBAND AD/NEO TRANSDUCER 1/50       |
| OXICLIQ-PI       | PEDIATRIC OXYGEN TRANSDUCER                  |
| OXI-P/I          | NELLCOR OXIBAND PED/INF TRANSDUCER 1/50      |

**Table 6. Radio Components**

| <b>Component Number</b> | <b>Object Description</b>                                     |
|-------------------------|---------------------------------------------------------------|
| 708920                  | Spot Vital Signs LXi Wireless Radio (A/B/G)                   |
| 703380                  | Spot Vital Signs LXi Wireless Radio Mounting Bracket          |
| 704622                  | Spot Vital Signs LXi Wireless Radio Mounting Instructions     |
| 4500-915                | Image Team 4600 2D Barcode Scanner Kit                        |
| 703878                  | Null Modem 9-Pin M/F 4-40                                     |
| 4500-920                | Kit, Spot Vital Signs LXi 802.11b Radio                       |
| 4500-921                | 802.11 b radio Directions for Use                             |
| 4500-922                | Kit, Wireless Spot Vital Signs LXi 802.11a/b/g Radio          |
| 4500-923                | 802.11a/b/g Radio Directions for Use                          |
| 4500-905                | Repair Tool CD, Spot Vital Signs LXi                          |
| 4500-906                | Spot Vital Signs LXi Configuration CD                         |
| 4500-907                | Spot Vital Signs LXi Firmware Loader CD                       |
| 4500-925                | Cable for wired connectivity                                  |
| 4500-926                | Cable for wired connectivity, Keyspan (USB to serial adapter) |
| 4500-927                | USB 2.0/serial cable kit                                      |
| 404684                  | Spot Vital Signs LXi Firmware Upgrade CD in case              |
| 404685                  | Spot Vital Signs LXi Radio Configuration CD in case           |



# 11

## Disassembly and Repair

This chapter provides the instructions for removing and replacing serviceable modules in the Spot Vital Signs LXi.

In general, re-assembly procedures are the reverse order of the disassembly procedures.



**WARNING** Electric shock hazard. There are no user-serviceable parts inside Spot Vital Signs LXi other than battery replacement (see “Battery Replacement” on page 46). Only an individual specifically trained and approved for the repair and/or verification of this specific Welch Allyn product may perform maintenance procedures specifically described in this manual. For service, refer the device to an Authorized Service Center.

**Note** Always disconnect the sealed lead-acid battery in the Spot Vital Signs LXi before performing any repair function.

### Tool list


| Description                                 | Part Number | Qty | Source            |
|---------------------------------------------|-------------|-----|-------------------|
| Blood Pressure Test Volume Repair Fixture   | 407672      | 1   | Welch Allyn       |
| Spot LXi Service Tool Software              | 4500-905    | 1   | Welch Allyn       |
| Service Manual, Spot Vital Signs LXi        | 704432      | 1   | Welch Allyn       |
| Spot Vital Signs LXi Directions for Use     | 705310      | 1   | Welch Allyn       |
| 802.11b Radio                               | 4500-920    | 1   | Welch Allyn       |
| 802.11b Radio Directions for Use            | 4500-921    | 1   | Welch Allyn       |
| 802.11a/b/g Radio                           | 4500-922    | 1   | Welch Allyn       |
| 802.11a/b/g Radio Directions for Use        | 4500-923    | 1   | Welch Allyn       |
| Tester, Calibration, 9600 Plus              | 01802-110   | 1   | Welch Allyn       |
| SureTemp Oral Probe                         | 02895-000   | 1   | Welch Allyn       |
| Blood Pressure Y-Tube, No Fittings 1/8 Tube | 5082-183    | 1   | Welch Allyn       |
| Welch Allyn Cuff                            | Soft-11-2MQ | 1   | Welch Allyn       |
| Welch Allyn Blood Pressure Hose, 5ft.       | 4500-30     | 1   | Welch Allyn       |
| Cable, USB A to 5 pin mini                  | 704889      | 1   | Welch Allyn       |
| DPAC Wireless Module                        | 706514      | 1   | Welch Allyn       |
| Needle Nose Pliers                          |             | 1   | Tool/Supply Store |
| Wire Cutter                                 |             | 1   | Tool/Supply Store |

| Description                                        | Part Number | Qty | Source                               |
|----------------------------------------------------|-------------|-----|--------------------------------------|
| Tweezers                                           |             | 1   | Tool/Supply Store                    |
| #2 Phillips Screwdriver                            |             | 1   | Tool/Supply Store                    |
| T10 Torx Driver                                    |             | 1   | Tool/Supply Store                    |
| 4" PVC Pipe                                        |             | 1   | Tool/Supply Store                    |
| Soldering Station                                  |             | 1   | Tool/Supply Store                    |
| Cable Tie Tool - T9921                             | GS2B        | 1   | Marsh Electronics<br>+1 800 877-8919 |
| *Setra Pressure Meter                              | 2270-01     | 1   | Setra +1 800 257 3872                |
| *Nitech Pressure Meter                             | 2000 in     | 1   | Nitech +1 800 547 6557               |
| Nellcor Sensor Ext Cable (4' Cable)                | Dec-4       | 1   | Nellcor +1 800 635 5267              |
| Nellcor Sensor Ext Cable (8' Cable)                | Dec-8       | 1   | Nellcor +1 800 635 5267              |
| Nellcor Sensor                                     | DS-100A     | 1   | Nellcor +1 800 635 5267              |
| Masimo SpO <sub>2</sub> Tester                     | 1795        | 1   | Masimo +1 800 326-4890               |
| Nellcor SRC-MAX SpO <sub>2</sub> Tester            | SRC-MAX     | 1   | Nellcor +1 800 635 5267              |
| Masimo Patient Cable                               | 1005        | 1   | Masimo +1 800 326-4890               |
| Masimo Adult SpO <sub>2</sub> Reusable Sensor      | 1269        | 1   | Masimo +1 800 326-4890               |
| Digital Volt Meter with 4 1/2 Digit Display        |             | 2   | Electronics Supply Store             |
| Power Supply: 0-20 Vdc adjustable with 0-3A output |             | 1   | Electronics Supply Store             |
| Wireless Router                                    |             | 1   | Electronics Supply Store             |
| Ethernet Cable                                     |             | 1   | Electronics Supply Store             |
| IBM compatible computer Windows 2000, XP, NT4      |             | 1   |                                      |
| * Only requires one of the pressure meters         |             |     |                                      |


### To disassemble Spot Vital Signs LXi:

1. Disconnect the power and all accessories from the Spot Vital Signs LXi.
2. Remove the four screws holding the battery door using a phillips-head screwdriver. Remove the battery door to expose the battery.







3. Tip the Spot Vital Signs LXi backward and slide the battery out. Disconnect the one-way connector.
4. Remove the two screws inside the battery housing that are identified with arrows molded into the housing and the two screws at the top of the back upper housing.
5. Hold the device together, lay the back housing on the bench, and carefully lift the top housing off.




6. Remove the screw and slightly lift the circuit board from the right facing side and disconnect the flex cable (J103).




7. Observe the locations of the two tubes before disconnecting. Lift and disconnect the wire harnesses and the pneumatic tubing. The two wire harnesses are a 6-wire harness that connects to J104 and a 2-wire harness that connects to J80 on the small, mounted blood pressure board.




8. Disconnect the red electrical connector and the gold ribbon connector. This is located on the main board. Lift up the tab on the white cable retainer so cable can slide out after it is unplugged.



9. Lift the flap on the beige connector (large ribbon cable) to remove the main board. This releases the serial board cable from connector J108.






**Note** This step only applies to devices configured with Masimo SpO<sub>2</sub> sensors. Cut jumper JP2 with diagonal cutting pliers twice so that it is open. Cut it on each end on the dotted lines as illustrated, so that it is an open connection.

## Key Pad Disassembly

1. Gently pull the button switch array out of the cover.
2. Align and push the new button switch array into place.




## SpO<sub>2</sub> Circuit Board Disassembly

**Note** To assure proper SpO<sub>2</sub> operation, replace the SpO<sub>2</sub> board using only the Welch Allyn specified part.


To assure patient electrical isolation, after the main board is nearly back in position, verify that the SpO<sub>2</sub> flex cable is freely floating in space and is not pressed up against the main board.

## Masimo


**Note** Masimo configurations ship with an EMI foil shield, copper tape, and Nomex insulator. These components are removed from the illustrations to show clarity of detail. These components must be reinstalled to comply with EMC requirements.




1. Lay the LCD flat on an ESD mat.
2. Find the double-stack of circuit boards and remove the four corners screws. Carefully lift the circuit board straight up.




3. Remove the three screws on the circuit board located behind the LCD.



4. Slightly rock the circuit board back and forth while lifting straight up.




5. Replace with a new circuit board.




## Nellcor

1. Lay the LCD flat on an ESD mat.
2. Remove the three screws on the circuit board located behind the LCD.



3. Slightly rock the circuit board back and forth while lifting straight up.





4. Replace with a new circuit board and set the dip switches to match the original board (dip 1 and 2: up / dip 3 and 4: down).

## Blood Pressure Circuit Board Disassembly

**Note** When replacing the NIBP board see “NIBP board Initialization” on page 62 for instructions.

Remove the four screws and carefully lift the circuit board straight up. There is one stationary connector.






## LCD Inverter Ballast Board Disassembly


This is a small board held in place with two-sided tape. It is located towards the bottom of the main board. Remove the connector at each end and pull straight up while taking care not to damage the board because the adhesive tape has considerable sticking strength.

## LCD Disassembly

1. Remove the main board. At the bottom of this board is a circuit board with a pink connector at the bottom.



2. Pull the brown tabs of the top connector forward and remove the ribbon cable. Inspect connector tabs for wear and replace if worn.
3. Turn the board over and pull back the soft rubber holders. Pull the cable from the lug to remove the grounding strap and remove the LCD.



4. Replace with a new LCD.


## Pump Disassembly

**Note** For proper blood pressure operation, replace the pump using only the Welch Allyn specified part.

To assure patient electrical isolation, route the pump wires through the rear housing clip feature (near the pump terminals) and held in place with a tiewrap to the holes in the rear housing at the back of the battery compartment.

To assure patient electrical isolation, do not modify the length of the pump wires.


1. Snip the wire ties (2 places).



2. Disconnect the pneumatic tubing and unhook the wiring.
3. Replace with a new pump.

## RS232-Communication Circuit Board Disassembly

1. Remove the pneumatic connector and associated hoses from the pump and the two screws in the rear housing. Remove the associated hoses.
2. Release the catch and remove the cover and inner handle assembly.
3. Remove the connection port panel.
4. Remove the two screws inside the back housing and the ribbon cable that hold the RS232 board in place. Take care since the connector is easily broken.
5. Slide the board straight out and remove the wire harness.



6. Replace with a new circuit board.

## Fan and Power Circuit Board Disassembly

**Note** To prevent buildup of hydrogen gas, replace the fan using only the Welch Allyn specified part.

This is a small board held in place with two-sided tape and one connector. Remove the screw before disconnecting the connector.

## Thermometry Circuit Board Disassembly



**Caution** To assure patient electrical isolation, trap the flex cable to the thermometer pod behind the clip on the main housing adjacent to the main board connection point.

**Caution** To assure patient electrical isolation, verify that the correct flex cable is used to connect the thermometer pod to the main board. The Braun and SureTemp Plus flex cables are slightly different in length and are not interchangeable.


**Caution** To protect the user from high voltage, properly install the thermometer pod insulating paper separating the LCD ballast wires from the housing crack.

## SureTemp Plus




**Caution** Do not attempt to repair or clean the solder joints on either SureTemp Plus board. Incorrect flux or technique can degrade thermometer accuracy.


1. Remove the two screws located at the thermometer housing with a T-10 torx wrench. Save the Nomex.




2. Locate the path of the notched tabs on the housing. Slide the housing up and out following this path and set the case aside.



3. Remove the four screws from the mounting plate and save the Nomex. This exposes the circuit board.



4. Remove the three Torx screws (T10) on the board (one on bottom and one on either side of the probe well cover) and lift the board out of the housing and orient the part as shown.




5. Disconnect all connectors and replace with a new board.

To replace the temperature connector, remove the one screw holding the temperature board and connector. Verify the Nomex is properly mounted upon reassembly.


When replacing the SpO<sub>2</sub>/blood pressure housing, verify the orange and white twisted wire is pushed as far into the handle housing as possible. Failure to do so can cause the wires to stretch when the SpO<sub>2</sub>/blood pressure housing is inserted.

## Braun ThermoScan PRO 4000


1. Remove the two screws located at the thermometer housing with a T-10 torx wrench. Save the Nomex.



2. Locate the path of the notched tabs on the housing. Slide the housing up and out following this path and set the case aside.
3. Remove the four screws that hold the metal shield in place.



4. Remove the springs and posts. Carefully rock the board back and forth while lifting straight up.



5. Disconnect all connectors and replace with a new board.
6. Reassembly is the reverse of disassembly.



# 12

## Repair Test Specifications

---

This appendix refers to the Spot Vital Signs LXi without pneumatics (tubing and cuff), temperature probe, SpO<sub>2</sub> probe, or main battery attached, unless otherwise noted.

Use the Repair Software for the performing the tests on the Spot Vital Signs LXi. The standard test voltage, unless otherwise stated is 6.5 Vdc (+0/-0.25 Vdc). Unless otherwise stated, all calibrated volumes are  $\pm 10$ cc of the stated volume.

### General Unit Test

#### A-D Noise Test

The A-D Noise Test is defined as the amount of noise on the Spot Vital Signs LXi A-D pressure channel; over a 1 second sample time while 0.0 mmHg is applied to the Device pressure port. The maximum limit is 0.050 mmHg.

#### Leak Test

The Leak Test is defined as the amount of pressure drop that is recorded over a 15-second interval with a 100cc cylinder attached to the Spot Vital Signs LXi pressure port and that volume having a stabilized pressure of 250 mmHg. The maximum limit is 6 mmHg.

#### Inflation Test

The Inflation Test is defined as the amount of time the Spot Vital Signs LXi pump can inflate a 250 cc cylinder from 5 mmHg to 210 mmHg. The maximum limit is 7 seconds.

#### Dump Test

The Dump Test is defined as the amount of time it takes Spot Vital Signs LXi to deflate a 500 cc cylinder from 260 mmHg to less than 15 mmHg. The maximum limit is 10 seconds.

#### Pneumatic Accuracy Test

Enter the Internal Configuration Mode to perform this test (see "Temperature Functional Check" on page 77).

Perform an Auto-zero before starting this test (see "Temperature Functional Check" on page 77). The Pneumatic Accuracy test is defined as the comparison of Spot Vital Signs LXi's pressure measurement and applied pressure at 0, 50, 150, and 250 mmHg. The

specification for each reading is within  $\pm 3$  mmHg of the target pressure except for the 0 mmHg reading which is within  $\pm 1.0$  mmHg.

## Valve Control Test

Perform the valve control test using the Welch Allyn Spot LXi Service Tool, see “Valve control test” on page 67.

## Voltage Calibration

Calibrate the Spot Vital Signs LXi battery voltage measurement circuit at  $5.5 + 0.1$  V DC.

After performing a successful battery measurement circuit calibration, Spot Vital Signs LXi's memory will store a calibration signature of up to four printable characters.

## Blank Mode Current Test

The Blank Mode Current Test is defined as the amount of current drawn through the battery terminals. All LCD segments, the back light, and the SpO<sub>2</sub> mode are all turned off. The limit is 200 mA maximum.

## Back Light (Idle) Current Test

The Back Light Current Test is defined as the amount of current drawn through the battery. All LCD segments and the back light are turned on while the SpO<sub>2</sub> mode is turned off. The limit is 400 mA maximum.

## Valve/Pump Mode Current Test

The Valve/Pump Mode Current Test is defined as the amount of current that is drawn through the battery terminals; Spot Vital Signs LXi is in the Blank Mode while the valve and pump are actuated to on. The limit is 2.0 A maximum.

## Interface Test

When given the proper commands, the Display Window will turn on or off all segments and turn on or off the back light. When given the proper command, the buzzer sounds to verify operation. The user determines the pass/fail criteria.

# Temperature Option Requirements

## Accuracy Test

Verify the accuracy of the temperature module is within  $\pm 0.2^\circ$  F for readings with a normal temperature of  $97.3^\circ$  F ( $36.3^\circ$  C) using a Cal Key (5200-25). Verify Spot Vital Signs LXi can read a temperatures of  $96.4^\circ$  F ( $35.8^\circ$  C) and  $106^\circ$  F ( $41.1^\circ$  C) within  $\pm 0.3^\circ$  F/ $\pm 0.2^\circ$  C using a Welch Allyn 9600 Plus Calibrator.



## Temperature Probe Test

Spot Vital Signs LXi displays “ORL” after you remove the blue thermometer probe from the blue probe well.

## SpO<sub>2</sub> Option Requirements

These tests are only valid on Spot Vital Signs LXi models with the Masimo or Nellcor SpO<sub>2</sub> option.

## SpO<sub>2</sub> Functional Test

See [“To begin functional verification of the Masimo SpO<sub>2</sub>”](#) or [“To begin functional verification of the Nellcor SpO<sub>2</sub>”](#) on page 87 for the functional test.

## SpO<sub>2</sub> Mode Current Test

The SpO<sub>2</sub> Mode Current Test is defined as the amount of current, less the Blank Mode Current, that is drawn through the battery terminals, placing Spot Vital Signs LXi in the Blank Mode, actuating the SpO<sub>2</sub> mode and applying any SpO<sub>2</sub> signal to the device. The maximum limit is 120 mA maximum.

## Fail Safe Test

### Over Pressure Test


Verify that Spot Vital Signs LXi can detect over pressure on the pneumatic system between 296 mmHg and 329 mmHg.

### Over 15 mmHg

Verify that Spot Vital Signs LXi can detect static pressure over 15 mmHg for 180 seconds.



Interconnect Diagram 1









Advancing Frontline Care™

Reorder No. 4500-89  
Material No. 718440 Ver. A